**Enzootic Bovine Leukosis Virus (BLV)
Antibody Test Kit**

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Enzootic Bovine Leukosis Virus (BLV) Antibody Test Kit

For veterinary use only.

Name and Intended Use

IDEXX Leukosis Serum Screening is IDEXX's enzyme immunoassay for the detection of Antibodies directed against Enzootic Bovine Leukosis Virus (BLV) in individual serum and plasma samples from bovines, and pool of sera from bovines (maximum 10).

General Information

Enzootic Bovine Leukosis is an infectious lymphoproliferative disease in cattle which occurs throughout the world. The disease is caused by an exogenous delta-retrovirus, Enzootic Bovine Leukosis Virus which established a persistent infection in a sub-population of B lymphocytes by integration of proviral DNA at a number of sites on the cellular DNA. The majority of infected cattle remain healthy for life although approximately 30% develop a persistent lymphocytosis and a small portion (up to 10%) develops lymphoid tumors. The disease affects mainly dairy herds and spreads predominantly through horizontal transmission by exposure to blood or secretions containing infected lymphocytes. As no treatment or vaccine is available, the eradication programs are based on identification and elimination of infected animals which is realised mainly by detection of anti-viral Antibodies. The infected cattle produce specific Antibodies to the major viral proteins from an early stage of infection. Initially, the agar-gel immunodiffusion (AGID) test was widely used. However, the sensitivity of AGID test is limited and instances have been reported where infected animals failed to produce a detectable Antibody response. For this reason, ELISA method, which is simple, rapid, and sensitive method, is used for surveillance of herds.

Descriptions and Principles

Microplates are coated with BLV Antigen. Samples to be tested are diluted and incubated in the wells. Upon incubation of the test sample in the coated wells, BLV specific Antibodies form immune-complexes with BLV-Antigens. After washing away unbound material, an anti-bovine Antibody enzyme Conjugate is added which binds to any Antigen-Antibody immune-complex. Unbound Conjugate is washed away and enzyme Substrate (TMB) is added. In presence of the enzyme, the Substrate is oxidized and develops a blue compound becoming yellow after blocking. Subsequent color development is directly related to the amount of Antibody to BLV present in the test sample. The result is obtained by comparing the sample Optical Density with the Positive Control mean Optical Density.

Reagents**Volume**

		5	10
1	BLV Antigen Coated Plate		
2	Positive Control	1 x 1.0 mL	1 x 1.0 mL
3	Negative Control	1 x 1.0 mL	1 x 1.0 mL
4a	Conjugate Concentrate (100X)	1 x 1.5 mL	1 x 1.5 mL
4b	Dilution Buffer N.1	1 x 120 mL	1 x 120 mL
5	Dilution Buffer N.2	1 x 120 mL	2 x 120 mL
A	TMB Substrate N.13	1 x 60 mL	1 x 120 mL
B	Stop Solution N.3	1 x 60 mL	1 x 120 mL
C	Wash Concentrate (10X) N.2	2 x 100 mL	4 x 100 mL

Note: See table at the end of the insert for a description of symbols used on the insert and labels of this kit.

Storage

Store the reagents at 2–8°C. Reagents are stable until expiration date, provided they have been stored properly.

Materials Required but Not Provided

- Precision pipettes or multiple-delivery pipetting devices suitable for delivering 10 to 1000 μ L
- Disposable pipette tips
- Graduated cylinder for wash solution
- 96-well Microplate reader (equipped with 450 nm filter)
- Microplate washer (manual, semi-automatic or automatic system)
- Use only distilled or deionized water for preparation of the reagents used in the test
- Vortex or equivalent
- Microplate covers (lid, aluminium foil or adhesive)
- Centrifuge (2000 x g)
- Humid Chamber / Incubator capable of maintaining a temperature of 37°C (\pm 3°C)
- Microplate shaker
- Uncoated microplates or tubes for sample preparation

Precautions and Warnings

- Handle all biological material as potentially infectious.
- Wear protective gloves / protective clothing / eye or face protection when handling samples and reagents.
- Refer to the product Material Safety Data Sheet for additional information.
- See the end of this insert for reagent hazard and precaution warnings.

Laboratory Practices

- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing, and washing throughout this procedure are necessary to maintain precision and accuracy. Use a separate pipette tip for each sample and control.
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plastic ware.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional, and national regulations.
- Care should be taken to prevent contamination of kit components. Do not pour unused reagents back into containers.
- Do not use kit past expiration date and do not intermix components from kits with different lot numbers.

Preparation of Reagents

Controls and Samples

The Dilution of the Samples and Controls depends on the selected incubation protocol and should be done in a pre-dilution plate or in tubes:

- Short protocol (1 hour (\pm 5 min.) at $+37^{\circ}\text{C}$ (\pm 3°C)): pre-dilute Controls and Samples 1:20 with Dilution Buffer N.2.
- Overnight protocol (16–24 hours at $2-8^{\circ}\text{C}$): pre-dilute Controls and Samples 1:50 with Dilution Buffer N.2.

Wash Solution

The Wash Concentrate (10X) N.2 must be diluted 1:10 with distilled/deionized water before use (e.g. 30 ml of Wash Concentrate (10X) N.2 in 270 ml of distilled water). This solution is hereafter called "Wash Solution".

Note: The Wash Concentrate (10X) N.2 should be brought to $18-26^{\circ}\text{C}$ and well mixed to ensure dissolution of any precipitated salts. The Wash Solution is stable for up to 3 days when stored at $2-8^{\circ}\text{C}$.

Conjugate

The Conjugate Concentrate (100X) must be diluted 1:100 in the Dilution Buffer N.1.

Note: the diluted Conjugate solution is stable for up to 8 hours at $18-26^{\circ}\text{C}$.

Test Procedure

All reagents must be allowed to come to 18–26°C before use. Mix reagents by gentle inverting or swirling.

1 Obtain coated microplates and record the position of each sample.

2 Distribute samples and controls:

Short protocol: 1 hour (± 5 min.) at +37°C (± 3 °C):

- Dispense 100 μL of DILUTED Negative Control into one appropriate well.
- Dispense 100 μL of DILUTED Positive Control into two appropriate wells.
- Dispense 100 μL of DILUTED samples into remaining wells.
- Homogenize contents of the wells using a microplate shaker.
- Cover the microplate and incubate for 1 hour (± 5 min.) at +37°C (± 3 °C).

Overnight protocol : 16–24 hours at 2–8°C:

- Dispense 100 μL of DILUTED Negative Control into one appropriate well.
 - Dispense 100 μL of DILUTED Positive Control into two appropriate wells.
 - Dispense 100 μL of DILUTED samples into remaining wells.
 - Homogenize contents of the wells using a microplate shaker.
 - Cover the microplate and incubate for 16–24 hours at 2–8°C.
-

3 Remove the solution and wash each well with approximately 300 μL of Wash Solution three times. Avoid plate drying between plate washings and prior to the addition of the next reagent. Tap each plate onto absorbent material after the final wash to remove any residual wash fluid.

4 Dispense 100 μL of DILUTED Conjugate into each well.

5 Cover the microplate and incubate for 30 minutes (± 3 min.) at +37°C (± 3 °C).

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8 Incubate 20 minutes (± 3 min.) at 18–26°C away from direct light.

9 Dispense 100 μL of Stop Solution N.3 into each well.

10 Measure and record Optical Densities values of samples and Controls at 450 nm.

11 Calculate results.

Note: When using robotics, incubation of microplates in an incubation chamber allows working without plate covers. Use of robots is also not compatible with gentle microplate tapping or wiping. Plates can be held up to 1 hour in the dark prior to reading.

12 Calculations:

Controls

$$PC\bar{x} = \frac{PC1 (A450) + PC2 (A450)}{2}$$

Validity criteria

$$PC\bar{x} \geq 0.350$$

$$PC\bar{x} : NC(A450) \geq 3.00$$

For invalid assays, technique may be suspect and the assay should be repeated following a thorough review of the package insert.

Samples

$$S/P \% = \frac{\text{Sample (A450)} - NC (A450)}{PC\bar{x} - NC(A450)}$$

13 Interpretation:

Individual Sera and Plasma samples

Negative

$$S/P \% \leq 60$$

Positive

$$S/P \% > 60$$

Pools of Sera

Negative

$$S/P \% \leq 40$$

Positive

$$S/P \% > 40$$

Note: IDEXX has instrument and software systems available which calculate results and provide data summaries.

For technical assistance:

IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4895

IDEXX Europe Tel: +800 727 43399

Contact your IDEXX area manager or distributor or visit our website: idexx.com/contactlpd

IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries

© 2014 IDEXX Laboratories, Inc. All rights reserved.

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WARNING / ATTENTION / ATENCIÓN / ATENÇÃO / ACHTUNG

H316 / H319 / P280 / P332+P313 / P337+P313

Positive control / Negative Control – Causes mild skin irritation. Causes serious eye irritation. Wear protective gloves/eye protection/face protection. If skin irritation occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention.

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H319 / P280 / P337+P313 / EUH208

Dilution Buffer N.2 – Causes serious eye irritation. Wear eye protection/face protection. If eye irritation persists: Get medical advice/attention. Contains MIT. May produce an allergic reaction.

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H315 / H319 / P280 / P332+P313 / P337+P313

TMB Substrate – Causes skin irritation. Causes serious eye irritation. Wear protective gloves/eye protection/face protection. If skin irritation occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention.

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H302 / H315 / H317 / H319 / H335 / P280 / P333+P313 / P337+P313 / P363

Stop solution – Harmful if swallowed. Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. May cause respiratory irritation. Wear protective gloves/eye protection/face protection. If skin irritation or rash occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention. Wash contaminated clothing before reuse.

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IDEXX Montpellier SAS
IDEXX Europe B.V.