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06-02212-05

IDEXX

Equine Infectious Anemia Virus Antibody Test Kit

For veterinary use only.

Introduction

The purified antigen from Equine Infectious Anemia Virus (EIAV) can be used in the agar gel immunodiffusion (AGID) test to detect the presence of antibody against EIAV in the serum of infected horses. Nonspecific precipitin lines are virtually eliminated from the AGID test by using highly purified EIAV antigens. These non-specific precipitin lines can interfere with the reading of the test.

Description

The movement of antigen and antibody toward one another in an agar gel can result in the formation of precipitin lines. By applying this principle, the AGID test has been shown to reliably detect specific antibody formed after 10 to 30 days of infection with EIAV. Precipitin lines of identity can be visualized in strong positive sera as continuing from the lines between the control serum and center antigen wells. Weakly positive sera are indicated as a bend of the control serum lines towards the antigen well. No precipitin lines of identity form if the serum is negative.

Reagents (Contains reagents for up to 230 Tests)

Volume

		Volume
1	Equine Infectious Anemia Antigen, preserved with sodium azide (Black Cap).	1 x 3.9 mL
2	Equine Infectious Anemia Positive Control Serum, preserved with sodium azide (Red Cap).	1 x 11.7 mL

NOTE: See table at the end of the insert for a description of international symbols used on the labels of this kit.

Precautions and Warnings

1. Handle all samples and reagents as if capable of transmitting EIAV.
2. Dispose of contents in accordance with local, regional, and national regulations.
3. Do not use kit past expiration date and do not intermix components from kits with different serial numbers.
4. In vitro use only.

Storage


Store the reagents at 2–8°C. Reagents are stable until expiration date, provided they have been stored properly.

Specimen Collection and Preparation for Testing

Use only horse serum for test specimens. Specimens may be stored at 2–8°C up to five days. If longer storage is desired, store at -20°C. The presence of turbidity, hemolysis or visible indications of bacterial growth may interfere with the performance and accuracy of the test. Do not test specimens that exhibit bacterial growth.

Plate Preparation

Prepare AGID plates according to the protocols of national and regional animal health authorities. Refer to OIE Manual of Diagnostic Tests and Vaccines for Terrestrial Animals 2012, Chapter 2.5.6 – www.oie.int

Serum and Antigen Placement in Wells (Pattern I)

1. Completely fill the three alternate outside wells (see wells labeled 1, 2 and 3) with each of the three test sera. Avoid any overflow into the agar surface. Alternative patterns may be used where allowed by national or regional animal health authorities.
2. Fill the center well with 50 μL of purified antigen in the same manner (see well “A”).
3. Fill the three remaining outside wells with 50 μL of Positive Control Serum in the same manner (see wells B).
4. Incubate plates for 24–48 hours at 18–26°C in a humid chamber. Weak positive reaction lines may continue to develop for up to 48 hours.

Results

Pattern II shows different precipitin patterns in the AGID test after incubation:

1. Negative serum. Control precipitin lines (B wells) run straight into the negative sample well or bend slightly back towards the B wells (see well 1).
2. Weakly positive serum. Control precipitin lines (B wells) run into the sample well but bend away from the B wells and towards each other (see well 3).
3. Strongly positive serum. A precipitin line forms between the sample well and center well A which is continuous with the control precipitin lines between wells B and A (see well 2).
4. Very strongly positive serum. Control precipitin lines (B wells) bend towards each other before they reach the test serum wells. They may be connected by a broad diffuse precipitin line which will be close to the antigen (A) well. A sharper line can be observed if the serum is two-fold serially diluted and retested.
5. Weak positives may result from:
 - a. Serum taken during the incubation period of EIA. Retesting after two or three weeks may show stronger reaction.
 - b. If the horse continues to show no clinical signs and an unchanged weak positive AGID reaction, then the horse may be considered to be a carrier.
 - c. Weak positives may also be characteristic of nursing foals of infected mares. Retest of the foal after six months should show negative AGID results if the foal is uninfected. If the foal's mare (dam) is negative in the AGID test then the foal is probably infected.

Quality Control

If the positive control included in the kit does not react, do not use the kit. Call IDEXX Technical Services.

IDEXX Technical Services:

IDEXX USA Tel: 1 800 548 9997 or 1 207 556 4890

Fax: 1 800 328 5461 or 1 207 556 4826

IDEXX Europe Tel: 00800 727 43399

Fax: 00800 433 99329

U.S. Vet. License No. 313
Product Code 5515.22

*IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.
©2013 IDEXX Laboratories, Inc. All rights reserved.

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© 2013 IDEXX Laboratories, Inc. Alle Rechte vorbehalten.
IDEXX Europe B.V.