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BVDV/MD/BDV p80 Protein Antibody Test Kit

For veterinary use only

Name and Intended Use

IDEXX BVDV p80 Ab is IDEXX's enzyme immunoassay for the detection of Antibodies directed against p80 protein for diagnostic of Bovine Viral Diarrhea Virus (BVDV) and Mucosal Disease (MD) in individual serum, plasma and milk samples and in pools of serum (maximum 10) and tank milk samples from bovine origin and for diagnostic of Border Disease (BD) in individual serum and plasma samples and pools of serum samples (maximum 5) from sheep.

General Information

The Bovine Viral Diarrhoea Virus (BVDV) is a pestivirus, which causes two diseases: Bovine Viral Diarrhoea (BVD) and Mucosal Disease (MD). Bovine Viral Diarrhoea is a pathology induced by one of the two strains of the virus (cytopathogenic and non-cytopathogenic). The acute form, characterised by fever and diarrhoea, is transient, with high morbidity rates and low mortality rates. Adult animals can also be infected with an asymptomatic subclinical form. Mucosal Disease has low morbidity rates (1%), but high mortality rates. Most often it is characterized by ulcers at different levels of the digestive tract and diarrhoea that is often hemorrhagic. Many pathologies are associated with or aggravated by the BVDV, including, among others, respiratory diseases, slow development, congenital defects, etc... Mucosal Disease occurs in calves infected during gestation (Immunotolerant Persistently Infected (I.P.I.) animals). These animals were infected with a non-cytopathogenic strain by the transplacental way between the 42nd and the 120th day of gestation. This corresponds to a period when immunocompetence is being established in the foetus: foreign antigens present at this time are considered as self-antigens and no immune response is developed against them. Thus, the persistently infected animals do not produce antibodies against the strain they are infected with. Mucosal Disease is induced by mutation of the non cytopathogenic strain to a cytopathogenic strain. The main sources of infection are I.P.I. animals, which continuously produce and shed the virus, and, in a transient manner (during 10 days), animals recently contaminated with a primary BVDV infection. Transmission of the virus can be oral nasal, conjunctival, genital or transplacental. The presence of BVDV in a herd can be detected with serological screening, which reveals the presence of animals with specific antibodies. However, it does not enable persistently infected animals to be detected. Previously based on seroneutralisation methods, this new serological screening test uses the ELISA method. This test can be rapidly implemented. It is reliable and particularly well suited to the analysis of a large number of samples and can be used for tank milk analysis. It is based on the principle of competition between the bovine antibodies and a Peroxidase coupled monoclonal anti-p80-antibody "WB112".

This kit is designed for the detection of Antibodies directed against p80 protein for diagnostic of Bovine Viral Diarrhea Virus (BVDV) and Mucosal Disease (MD) in individual serum plasma and milk samples and in pools of serum (maximum 10) and tank milk samples from bovine origin, and for diagnostic of Border Disease (BD) in individual serum and plasma samples and pools of serum samples (maximum 5) from sheep.

Descriptions and Principles

Microplates are coated with p80 proteins attached to the wells via a specific WB103 antibody. Samples to be tested are diluted and incubated in the wells. Upon incubation of the test sample in the coated wells, p80 protein specific Antibodies form immune-complexes with the p80 protein. After washing away unbound material, an anti-p80 protein Antibody enzyme Conjugate is added. In presence of the p80 protein-Antibody immune-complex, the Conjugate is blocked from binding to its corresponding epitops on the microplate. Conversely, in the absence of p80 protein-Antibodies in the test sample, the Conjugate is free to bind to its corresponding epitops on the microplate. Unbound Conjugate is washed away and an enzyme Substrate (TMB) is added. In presence of the enzyme, the Substrate is oxidized and develops a blue compound becoming yellow after blocking. Subsequent color development is inversely proportional to the amount of anti-p80 protein Antibodies in the test sample (see “Calculations” and “Interpretation of results”).

The diagnostic relevance of the result is obtained by comparing the sample Optical Density with the Negative Control mean Optical Density (see “Calculations” and “Interpretation of results”).

Reagents

Store all reagents at 2–8°C.

Reagents	Volume	
1	BVDV/MD/BDV p80 Protein Coated Plates	5
2	Positive Control	2 mL
3	Negative Control	2 mL
4a	(Anti-p80 HRPO) Conjugate Concentrate	0.75 mL
4b	Dilution Buffer N.1	120 mL
5	Dilution Buffer N.9	120 mL
A	TMB Substrate N.9	60 mL
B	Stop Solution N.3	60 mL
C	Wash concentrate (20X)	100 mL

NOTE: see table on page 45 for the description of international symbols used on the kit labels.

Materials Required but Not Provided

- Centrifuge (capacity 2000 x g)
- Vortex or equivalent
- Precision Micropipettes and Multi-dispensing micropipettes (reagents volumes listed in the “Test Protocol” require pipette precision less than or equal to 5%)
- Disposable pipette tips
- Microplate shaker
- Distilled water or deionized water
- Microplate washer (manual, semi-automatic or automatic system)
- Microplate covers (lid, aluminium foil or adhesive)
- 96-well Microplate reader equipped with 450 nm filter

Precautions and Warnings for Users

- Do not pipette by mouth.
- Wear protective gloves/protective clothing/eye protection/face protection.
- Controls, TMB-Substrate and Wash concentrate solutions can cause eye irritation.
- Stop Solution can cause severe skin burns and eye damage.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional and national regulations.

Preparation of Reagents

Wash Solution

The Wash Concentrate (20X) must be diluted 1:20 with distilled/deionized water before use (e.g. 15 ml of Wash Concentrate (20X) in 285 ml of distilled water). This solution is hereafter called "Wash Solution".

Note: The Wash Concentrate (20X) should be brought to 18–26°C and well mixed to ensure dissolution of any precipitated salts. Wash Solution is stable for up to 3 days when stored at 2–8°C.

Conjugate

The Conjugate Concentrate must be diluted 1: 100 in the Dilution Buffer N.1.

Note: Diluted Conjugate solution is stable for up to 8 hours at 18–26°C.

Milk samples

Defatted as well as whole milk samples can be tested.

Test Procedure

All reagents must be allowed to come to 18–26°C before use.

Reagents should be mixed by gentle swirling or vortexing. Use a separate pipette tip for each sample. Controls may be dispensed anywhere on the microplate.

Obtain coated microplates and record the position of each sample on a worksheet.

1. Dispense Dilution Buffer N. 9 Controls and samples:

- a) Bovine Individual Serum and Plasma samples: Short incubation [1 hour (\pm 5 min.) at 18–26°C]
 - Dispense 50 μ L of Dilution Buffer N.9 in each well.
 - Dispense 50 μ L of Negative Control into two appropriate wells.
 - Dispense 50 μ L of Positive Control into one appropriate well.
 - Dispense 50 μ L of sample into remaining wells (1 well per sample).
 - Homogenize contents of the wells using a microplate shaker.
 - Cover the microplate (with a lid, aluminium foil or adhesive plate cover) and incubate for 1 hour (\pm 5 min.) at 18–26°C.
- b) Bovine Individual Serum and Plasma samples: Overnight incubation [16–24 hours at 2–8°C]
 - Dispense 90 μ L of Dilution Buffer N.9 in each well.
 - Dispense 10 μ L of Negative Control into two appropriate wells.
 - Dispense 10 μ L of Positive Control into one appropriate well.
 - Dispense 10 μ L of sample into remaining wells (1 well per sample)
 - Homogenize contents of the wells using a microplate shaker.
 - Cover the microplate (with a lid, aluminium foil or adhesive plate cover) and incubate for 16–24 hours at 2–8°C.

- c) Bovine Pool Serum samples (Maximum 10): Overnight incubation [16–24 hours at 2–8°C]
- Dispense 50 μL of Dilution Buffer N.9 in each well.
 - Dispense 50 μL of Negative Control into two appropriate wells.
 - Dispense 50 μL of Positive Control into one appropriate well.
 - Dispense 50 μL of sample into remaining wells (1 well per sample)
 - Homogenize contents of the wells using a microplate shaker.
 - Cover the microplate (with a lid, aluminium foil or adhesive plate cover) and incubate for 16–24 hours at 2–8°C.
- Note: Analysis of pools of serum samples from bovine origin is not approved for use in Germany.
- d) Bovine Individual and Tank Milk samples: Short incubation [2 hours (\pm 5 min.) at 18–26°C]
- Dispense 50 μL of Dilution Buffer N.9 in three appropriate wells.
 - Dispense 50 μL of Negative Control into two appropriate wells (containing 50 μL of Dilution Buffer N.9).
 - Dispense 50 μL of Positive Control into one appropriate well (containing 50 μL of Dilution Buffer N.9).
 - Dispense 100 μL of sample into remaining wells (1 well per sample).
 - Homogenize contents of the wells using a microplate shaker.
 - Cover the microplate (with a lid, aluminium foil or adhesive plate cover) and incubate for 2 hours (\pm 5 min) at 18–26°C.
- Note: Tank milk samples analysis is not approved for use in Germany.
- e) Bovine Individual and Tank Milk samples : Overnight incubation [16–24 hours at 18–26°C]
- Dispense 100 μL of Dilution Buffer N.9 in each well.
 - Dispense 100 μL of Negative Control into two appropriate wells.
 - Dispense 100 μL of Positive Control into one appropriate well.
 - Dispense 100 μL of sample into remaining wells (1 well per sample).
 - Homogenize contents of the wells using a microplate shaker.
 - Cover the microplate (with a lid, aluminium foil or adhesive plate cover) and incubate for 16–24 hours at 18–26°C.
- Note: Tank milk samples analysis is not approved for use in Germany.
- f) Sheep Individual Serum and Plasma samples and Pool of Serum samples (Maximum 5): Overnight incubation [16–24 hours at 2–8°C]
- Dispense 50 μL of Dilution Buffer N.9 in three appropriate wells.
 - Dispense 75 μL of Dilution Buffer N. 9 in remaining wells.
 - Dispense 50 μL of Negative Control into two appropriate wells (containing 50 μL of Dilution Buffer N. 9).
 - Dispense 50 μL of Positive Control into one appropriate well (containing 50 μL of Dilution Buffer N. 9).
 - Dispense 25 μL of sample per well (containing 75 μL of Dilution Buffer N. 9).
 - Homogenize contents of the wells using a microplate shaker.
 - Cover the microplate (with a lid, aluminium foil or adhesive plate cover) and incubate for 16–24 hours at 2–8°C.
- Note: Analysis of pools of serum samples from ovine origin is not approved for use in Germany.

2. Wash each well with approximately 300 μL of Wash Solution 3 to 5 times. Aspirate the liquid contents of all well after each wash. Following the final aspiration, firmly tap residual wash fluid from each microplate onto absorbent material. Avoid microplate drying between washes and prior to the addition of the next reagent.
- Note: If skimmed milk (or milk taken from under the cream layer) is used, this type of washing is sufficient. With full fat milk, it is recommended to modify this method of washing by addition of a soak step for 1 minute per wash cycle. This facilitates the elimination of fat particles that are likely to fix in a non-specific way the conjugate in the next step. Thorough washing is essential for optimal results.
3. Dispense 100 μL of diluted Conjugate into each well.
 4. Cover the microplate (with a lid, aluminium foil or adhesive plate cover) and incubate for 30 minutes (\pm 3 min.) at 18–26°C.
 5. Wash each well with approximately 300 μL of Wash Solution 3 times. Aspirate the liquid contents of all well after each wash. Following the final aspiration, firmly tap residual wash fluid from each microplate onto absorbent material. Avoid microplate drying between washes and prior to the addition of the next reagent.
 6. Dispense 100 μL of TMB Substrate N.9 into each well.
 7. Incubate 20 minutes (\pm 3 min.) at 18–26°C in a dark place.
 8. Dispense 100 μL of Stop Solution N.3 into each well. Shake the microplate by gentle tapping. Wipe carefully the underside of the microplate.
 9. Blank the microplate reader on air.
 10. Measure and record Optical Densities values of samples and Controls at 450 nm.
 11. Calculate results.

Note: When using robotics, incubation of microplates in an incubation chamber allows working without plate covers. Use of robots is also not compatible with gentle microplate tapping or wiping. Plates can be held up to 1 hour in the dark prior to reading.

Results

For the assay to be valid, the Negative Control mean ($\text{NC}\bar{x}$) must be greater than or equal to 0.800. In addition, the S/N percentage of the Positive Control should be less than 20%.

NOTE: IDEXX has instrument and software systems available which calculate means and % values and provide data summaries.

Calculation

Negative Control Mean

$$NC\bar{x} = \frac{NC1 A_{450} + NC2 A_{450}}{2}$$

Calculation for test samples

$$S/N (\%) = 100 \times \frac{\text{Sample } A_{450}}{NC\bar{x}}$$

Interpretation of Results

Serum samples

BVD/MD diagnostic for bovine Individual Serum and Plasma samples and BD diagnostic for individual serum and plasma and pooled serum samples from sheep:

- Samples with S/N % greater than or equal to 50 % are considered Negative.
- Samples with S/N % greater than 40 % and less than 50 % are considered doubtful and must be retested.
- Samples with S/N % less than or equal to 40 % should be considered Positive.

BVD/MD diagnostic for bovine Pool of Serum samples :

- Samples with S/N % greater than or equal to 60 % are considered Negative.
- Samples with S/N % greater than 50 % and less than 60 % are considered doubtful and must be retested.
- Samples with S/N % less than or equal to 50 % should be considered Positive.

Notes:

- Positive pools should be confirmed by analyzing the individual samples composing the pool.
- Analysis of pools of serum samples is not applicable in Germany.

Milk samples

Qualitative interpretation - BVD/MD diagnostic for bovine Individual and Tank Milk samples:

- Individual and Tank Milk Samples with S/N % greater than or equal to 80 % are considered as samples coming from a negative animal or herd.
- Individual and Tank Milk Samples with S/N % less than 80 % should be considered as samples coming from an animal or herd which has been in contact with BVDV/MD.

Semi-quantitative interpretation - BVD/MD diagnostic for bovine Tank Milk samples:

- Samples with S/N % greater than or equal to 80 % are considered as coming from a herd with less than 10 % BVD/MD within herd sero-prevalence.
- Samples with S/N % greater than 45 % and less than 80 % are considered as coming from a herd with a 10 % to 30 % BVD/MD within herd sero-prevalence.
- Samples with S/N % less than or equal to 45 % are considered from herd with BVD/MD prevalence greater than 30 %.

Note:

- Analysis of Tank milk samples is not applicable in Germany.

👁 = *Modification in the using instructions*

Summarized Test Procedure

IDEXX strongly recommends that you read the complete instructions carefully before using the test the first time.

Steps	Action
1. Preparation of reagents	Wash Concentrate (20X) must be diluted 1:20 with distilled/deionized water before use Concentrated Conjugate must be diluted 1: 100 in the Dilution Buffer N. 1
2. Distribution and incubation of samples	<p>1. Dispense Dilution Buffer N. 9 Controls and samples:</p> <p><u>Bovine Individual Serum and Plasma samples: Short incubation [1 hour (\pm 5 min.) at 18–26°C]</u></p> <ul style="list-style-type: none"> • Dispense 50 μL of Dilution Buffer N.9 in each well. • Dispense 50 μL of Negative Control into two appropriate wells. • Dispense 50 μL of Positive Control into one appropriate well. • Dispense 50 μL of sample into remaining wells (1 well per sample). • Homogenize contents of the wells using a microplate shaker. • Cover the microplate and incubate for 1 hour (\pm 5 min.) at 18–26°C. <p><u>Bovine Individual Serum and Plasma samples: Overnight incubation [16–24 hours at 2–8°C]</u></p> <ul style="list-style-type: none"> • Dispense 90 μL of Dilution Buffer N.9 in each well. • Dispense 10 μL of Negative Control into two appropriate wells. • Dispense 10 μL of Positive Control into one appropriate well. • Dispense 10 μL of sample into remaining wells (1 well per sample) • Homogenize contents of the wells using a microplate shaker. • Cover the microplate and incubate for 16–24 hours at 2–8°C. <p><u>Bovine Pool Serum samples (Maximum 10): Overnight incubation [16–24 hours at 2–8°C]</u></p> <ul style="list-style-type: none"> • Dispense 50 μL of Dilution Buffer N.9 in each well. • Dispense 50 μL of Negative Control into two appropriate wells. • Dispense 50 μL of Positive Control into one appropriate well. • Dispense 50 μL of sample into remaining wells (1 well per sample) • Homogenize contents of the wells using a microplate shaker. • Cover the microplate and incubate for 16–24 hours at 2–8°C. <p><u>Bovine Individual and Tank Milk samples: Short incubation [2 hours (\pm 5 min.) at 18–26°C]</u></p> <ul style="list-style-type: none"> • Dispense 50 μL of Dilution Buffer N.9 in three appropriate wells. • Dispense 50 μL of Negative Control into two appropriate wells (containing 50 μL of Dilution Buffer N.9). • Dispense 50 μL of Positive Control into one appropriate well (containing 50 μL of Dilution Buffer N.9). • Dispense 100 μL of sample into remaining wells (1 well per sample). • Homogenize contents of the wells using a microplate shaker. • Cover the microplate and incubate for 2 hours (\pm 5 min) at 18–26°C. <p><u>Bovine Individual and Tank Milk samples : Overnight incubation [16–24 hours at 18–26°C]</u></p> <ul style="list-style-type: none"> • Dispense 100 μL of Dilution Buffer N.9 in each well. • Dispense 100 μL of Negative Control into two appropriate wells. • Dispense 100 μL of Positive Control into one appropriate well. • Dispense 100 μL of sample into remaining wells (1 well per sample). • Homogenize contents of the wells using a microplate shaker. • Cover the microplate and incubate for 16–24 hours at 18–26°C. <p><u>Sheep Individual Serum and Plasma samples and Pool of Serum samples (Maximum 5): Overnight incubation [16–24 hours at 2–8°C]</u></p> <ul style="list-style-type: none"> • Dispense 50 μL of Dilution Buffer N.9 in three appropriate wells. • Dispense 75 μL of Dilution Buffer N. 9 in remaining wells. • Dispense 50 μL of Negative Control into two appropriate wells (containing 50 μL of Dilution Buffer N. 9). • Dispense 50 μL of Positive Control into one appropriate well (containing 50 μL of Dilution Buffer N.9). • Dispense 25 μL of sample per well (containing 75 μL of Dilution Buffer N. 9) • Homogenize contents of the wells using a microplate shaker. • Cover the microplate (with a lid, aluminium foil or adhesive plate cover) and incubate for 16–24 hours at 2–8°C.

3. Washing the plate	Wash each well with approximately 300 μ L of Wash Solution 3 to 5 times.
4. Conjugate distribution	Dispense 100 μ L of diluted Conjugate into each well.
5. Conjugate incubation	Cover the microplate and incubate for 30 minutes (\pm 3 min.) at 18–26°C.
6. Washing the plate	Wash each well with approximately 300 μ L of Wash Solution 3 times.
7. Substrate distribution	Dispense 100 μ L of TMB Substrate N.9 into each well.
8. Substrate incubation	Incubate 20 minutes (\pm 3 min.) at 18–26°C in a dark place.
9. Stopping the reaction	Dispense 100 μ L of Stop Solution N.3 into each well.
10. Measure the plate	Blank the microplate reader on air. Measure and record Optical Densities values of samples and Controls at 450 nm. Calculate results.
11. Interpretation (S/N percentage)	<p>Serum samples</p> <p>BVD/MD diagnostic for bovine Individual Serum and Plasma samples and BD diagnostic for individual serum and plasma and pooled serum samples from sheep:</p> <ul style="list-style-type: none">• Samples with S/N % greater than or equal to 50 % are considered Negative.• Samples with S/N % greater than 40 % and less than 50 % are considered doubtful and must be retested.• Samples with S/N % less than or equal to 40 % should be considered Positive. <p>BVD/MD diagnostic for bovine Pool of Serum samples :</p> <ul style="list-style-type: none">• Samples with S/N % greater than or equal to 60 % are considered Negative.• Samples with S/N % greater than 50 % and less than 60 % are considered doubtful and must be retested.• Samples with S/N % less than or equal to 50 % should be considered Positive. <p>Milk samples</p> <p>Qualitative interpretation - BVD/MD diagnostic for bovine Individual and Tank Milk samples:</p> <ul style="list-style-type: none">• Individual and Tank Milk Samples with S/N % greater than or equal to 80 % are considered as samples coming from a negative animal or herd.• Individual and Tank Milk Samples with S/N % less than 80 % should be considered as samples coming from an animal or herd which has been in contact with BVDV/MD. <p>Semi-quantitative interpretation - BVD/MD diagnostic for bovine Tank Milk samples:</p> <ul style="list-style-type: none">• Samples with S/N % greater than or equal to 80 % are considered as coming from a herd with less than 10 % BVD/MD within herd sero-prevalence.• Samples with S/N % greater than 45 % and less than 80 % are considered as coming from a herd with a 10 % to 30 % BVD/MD within herd sero-prevalence.• Samples with S/N % less than or equal to 45 % are considered from herd with BVD/MD prevalence greater than 30 %.

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For technical assistance:

Contact your IDEXX area manager or distributor
or visit: www.idexx.com/production/contact
IDEXX Technical Support: 00-800-727-43399

*IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

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IDEXX Montpellier SAS
IDEXX Europe B.V.