***Brucella abortus* Antibody Test Kit / Bovine Milk**

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For veterinary use only

Name and Intended Use

IDEXX Brucellosis Milk X2 ELISA Test Kit provides a rapid, simple, sensitive and specific method for detecting antibodies against *Brucella abortus* (*B. abortus*) in bovine individual and bulk milk samples.

- ☞ In Member States of the European Union, ELISA tests done on milk must be conducted on bulk milk samples as required in Annex C of EU Decision 2008/984/EC.

Descriptions and Principles

Microtiter plates are supplied precoated with inactivated antigen. Dilutions of the samples to be tested are incubated in the wells of these plates. Any antibody specific for *B. abortus* binds to the antigen in the wells and forms an antigen/antibody complex on the plate well surface. Unbound material is removed from the wells by washing. A peroxidase-labeled anti-ruminant IgG Conjugate is added, which binds to the ruminant antibodies complexed with the *B. abortus* antigen. Unbound Conjugate is removed by washing and the TMB Substrate is added to the wells. The degree of color that develops (optical density measured at 450 nm) is directly proportional to the amount of antibody specific for *B. abortus* present in the sample. The diagnostic relevance of the result is obtained by comparing the optical density (OD) that develops in wells containing the samples with the OD from wells containing the Positive Control.

Reagents		Volume
1	<i>B. abortus</i> Antigen Coated Plate	10
2	Positive Control (lyophilized)	1 x 5.0 mL
3	Negative Control (lyophilized)	1 x 5.0 mL
4	Conjugate	1 x 110 mL
A	TMB Substrate N.12	1 x 100 mL
B	Stop Solution N.3	1 x 100 mL
C	Wash Concentrate (10X)	1 x 480 mL
Other components: Zip lock bag		1

Note: See table at the end of the insert for a description of international symbols used on the labels of this kit.

Storage

Store the reagents at 2–8°C. Reagents are stable until expiration date, provided they have been stored properly.

Materials Required but Not Provided

- Precision micropipettes or multi-dispensing micropipettes
- Disposable pipette tips
- Graduated cylinder for wash solution
- 96-well microplate reader (equipped with 450 nm filter)
- Use only distilled or deionized water for preparation of the reagents used in the test.
- Microplate washer (manual, semi-automatic or automatic system)
- Microplate covers (lid, aluminium foil or adhesive)
- Humid Chamber/ Incubator capable of maintaining a temperature of +37°C ($\pm 3^{\circ}\text{C}$)
- Centrifuge (capacity 2000 x g)
- Microplate shaker

Precautions and Warnings for Users

- Handle all biological material as potentially infectious.
- The Substrate solution is irritating to eyes, respiratory system, and skin. Avoid contact with skin and eyes.
- Wear protective gloves/protective clothing/eye or face protection when handling samples and reagents
- Refer to the product Material Safety Data Sheet for additional information.

Laboratory Practices

- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing, and washing throughout this procedure are necessary to maintain precision and accuracy.
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plastic ware.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional, and national regulations.
- Care should be taken to prevent contamination of kit components. Do not pour unused reagents back into containers.
- Do not use kit past expiration date and do not intermix components from kits with different lot numbers.

Preparation of Reagents

Controls

Reconstitute Negative and Positive Controls with 5 ml of distilled/deionized water.

Notes: Reconstituted Controls must be aliquoted and kept frozen.

Wash Solution

Determine the amount of Wash Solution needed for washing the microtiter plates and diluting the samples and controls. Dilute the Wash Concentrate (10X) 1/10 with water (1 part Wash Concentrate with 9 parts water, e.g., 100 mL Wash Concentrate (10X) + 900 mL distilled water). When prepared under sterile conditions, the Wash Solution can be stored for one week at 2–8°C.

Preparation of Reagents

Milk samples must be skimmed milk samples (obtained after centrifugation).

Test Procedure

All reagents must be allowed to come to 18–26°C before use. Reagents should be mixed by gentle swirling or vortexing.

1. Dilute the Positive Control 1/4 in a separate tube with the Wash Solution.
2. Dispense 50 μL of Wash Solution into each well of the microtiter plate.
3. Add 50 μL of the DILUTED Positive Control into two appropriate wells.
Final dilution of positive control = 1/8
Add 50 μL of UNDILUTED Negative Control into two appropriate wells.
Add 50 μL of UNDILUTED skimmed milk samples into the appropriate wells.
Final dilution of samples and Negative Control milk = 1/2
4. Mix the contents within each well by gently shaking the microtiter plate briefly (a microtiter plate shaker might be used).
5. Cover the plate and incubate for:
 - Individual milk and bulk milk up to 125 cows: 60 minutes (± 5 min.) at +37°C (± 3 °C) or 14–18 hours at 2–8°C.
 - Bulk milk from 126 cows up to 250 cows: 14–18 hours at 2–8°C.With either options, plates should be hermetically sealed or covered and incubated in a humid chamber.
6. Empty or aspirate liquid content from the microwells and wash each well with approximately 300 μL of Wash Solution three times. Discard the liquid contents of all wells after each wash. Following the final fluid removal, firmly tap residual washfluid from each plate onto absorbent material. Avoid plate drying between washes and prior to the addition of the next reagent.
7. Dispense 100 μL of Conjugate into each well.
8. Cover the plate and incubate for 60 minutes (± 5 min.) at +37°C (± 3 °C). The plate(s) should be hermetically sealed or covered and incubated in a humid chamber.
9. Repeat Step 6.
10. Dispense 100 μL TMB Substrate N.12 into each well.
11. Incubate at 18–26°C for 15 minutes (± 1 min.).
12. Stop the color reaction by adding 100 μL of Stop Solution N.3 per well. The Stop Solution should be dispensed in the same order and at the same speed as the Substrate.
13. Read the results using a photometer at a wavelength of 450 nm.

Results

To validate the assay, the mean optical density (OD) of the Positive Control (PC \bar{x}) should not exceed 2.000 and the mean OD of the Negative Control (NC \bar{x}) should not exceed 0.500. The difference between the Positive and the Negative Control mean ODs (PC \bar{x} –NC \bar{x}) must be ≥ 0.300 . Make sure to read the plates within two hours after the addition of the Stop Solution.

Note: IDEXX has instrument and software systems available that calculate means and % values and provide data summaries.

Calculation

The ODs of samples tested in duplicates must be averaged. The mean OD of the Positive Control (PC \bar{x}) and the OD of the sample (Sample A450) are corrected by subtracting the mean OD of the Negative Control (NC \bar{x}):

Positive Control corrected OD:	Sample corrected OD:	S/P ratio for each sample:
PC \bar{x} – NC \bar{x}	Sample A450 – NC \bar{x}	S/P % = 100 x $\frac{\text{Sample A450} - \text{NC}\bar{x}}{\text{PC}\bar{x} - \text{NC}\bar{x}}$

Interpretation of Results

Interpretation for individual milks and bulk milk samples up to 125 cows:

S/P % (short incubation)	< 30 %	≥ 30 %
S/P % (long incubation)	< 60 %	≥ 60 %
Interpretation	Negative	Positive

Interpretation for bulk milk samples from 126 up to 250 cows:

S/P % (long incubation)	< 30 %	≥ 30 %
Interpretation	Negative	Positive

If a pool yields a positive result, it is recommended that test samples used to prepare the pool should be analyzed individually.

👁 = *Major change in the user instructions*

Summarized Test Procedure:

IDEXX strongly recommends that you read the complete instructions carefully before using the test for the first time.

Step	Action
1. Reagents	Reconstitute the lyophilized Control milks each with 5 ml of sterile water. Dilute the Wash Concentrate (10X) 1/10 with water to prepare the Wash Solution.
2. Sample Distribution	Dilute the Positive Control 1/4 in a separate tube with the Wash Solution. Dispense 50 μ L of Wash Solution into each well of the microtiter plate. Add 50 μ L of the DILUTED Positive Control into two appropriate wells. Final dilution of Positive Control = 1/8 Add 50 μ L of UNDILUTED Negative Control into two appropriate wells. Add 50 μ L of UNDILUTED skimmed milk samples into the appropriate wells Final dilution of samples and negative control milk = 1/2. Mix the contents within each well by gently shaking the microtiter plate briefly (a microtiter plate shaker might be used).
3. Sample Incubation	Cover the plate and incubate for: <ul style="list-style-type: none"> Individual milk and bulk milk up to 125 cows: 60 minutes (\pm5 min.) at +37°C (\pm3°C) or 14–18 hours at 2–8°C. Bulk milk from 126 cows up to 250 cows: 14–18 hours at 2–8°C in a humid chamber. With either options, plates should be hermetically sealed or covered and incubated in a humid chamber.
4. Washing the Plate	Empty or aspirate liquid content from the microwells and wash each well with approximately 300 μ L of Wash Solution three times.
5. Conjugate distribution	Dispense 100 μ L of Conjugate into each well.
6. Conjugate incubation	Cover the plate and incubate it for 60 minutes (\pm 5 min.) at +37°C (\pm 3°C). The plate(s) should be hermetically sealed or covered and incubated in a humid chamber.
7. Repeat step 4	
8. Substrate distribution	Dispense 100 μ L TMB Substrate N.12 into each well.
9. Substrate incubation	Incubate at 18–26°C for 15 minutes (\pm 1 min.).
10. Stopping the reaction	Stop the color reaction by adding 100 μ L Stop Solution N.3 per well.
11. Measure the plate	Read the results using a photometer at a wavelength of 450 nm.
12. Interpretation	Interpretation for individual milk and bulk milk samples up to 125 cows:
	S/P % (short incubation) < 30 % ≥ 30 %
	S/P % (long incubation) < 60 % ≥ 60 %
	Interpretation Negative Positive
	Interpretation for bulk milk samples from 126 up to 250 cows:
	S/P % (long incubation) < 30 % ≥ 30 %
Interpretation Negative Positive	
	If a pool yields a positive result, it is recommended that test samples used to prepare the pool should be analyzed individually.

For technical assistance:

Contact your IDEXX area manager or distributor or visit our website: idexx.com/contact/ptd
IDEXX Technical Support: +800 727 43399

*IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

©2013 IDEXX Laboratories, Inc. All rights reserved.

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©2013 IDEXX Laboratories, Inc. Alle Rechte vorbehalten.
IDEXX Switzerland AG
IDEXX Europe B.V.