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Avian Pneumovirus Antibody Test Kit

For veterinary use only

Name and Intended Use

IDEXX APV is IDEXX's enzyme immunoassay for the detection of antibody to avian pneumovirus (APV) in chicken and turkey serum.

General Information

Avian Pneumovirus (APV) is a pathogen of chickens and turkeys that can produce respiratory signs and reproductive failures as a primary pathogen, or as part of a multiagent syndrome known as Swollen Head Syndrome (SHS). Virus isolation of APV is a difficult task, so diagnosis through serology is one of the most important tools for disease control. This test is intended for quantitative monitoring chicken and turkey flocks that have been vaccinated with live and/or inactivated products containing A, B or C antigens, as well as naive (SPF) flocks for the detection of infected birds with any A, B or C serotype of APV.

Descriptions and Principles

This assay is designed to measure the relative level of antibody to APV in chicken and turkey serum. Viral antigen is coated on 96-well plates. Upon incubation of the test sample in the coated well, antibody specific to APV forms a complex with the coated viral antigens. After washing away unbound material from the wells, a conjugate is added that binds to any attached chicken and turkey antibody in the wells. Unbound conjugate is washed away and enzyme substrate is added. Subsequent color development is directly related to the amount of antibody to APV present in the test sample.

Reagents

Store all reagents at 2–8°C.

Reagents	Volume	
1	APV Antigen coated plates	5
2	Positive Control	1.9 ml
3	Negative Control	1.9 ml
4	(Anti-Chicken/Anti-Turkey HRPO) Conjugate	50 ml
5	Sample Diluent	235 ml
A	TMB Substrate	60 ml
B	Stop Solution	60 ml

NOTE: see table on page 28 for the description of international symbols used on the kit labels.

Materials Required but Not Provided

- Precision Micropipettes and Multi-dispensing micropipettes (reagents volumes listed in the “Test Protocol” require pipette precision of $\pm 5\%$)
- 96-well plate reader equipped with 650 nm filter
- Tubes for diluting samples
- Distilled or deionized water
- Device for the delivery and aspiration of wash solution

Precautions and Warnings for Users

- Handle all biological material as potentially infectious.
- Do not pipette by mouth.
- Do not eat, drink or smoke where specimens or kit reagents are being handled.
- The substrate solution is irritating to the eyes, respiratory system and skin. Avoid contact with skin and eyes.
- Wear protective gloves/protective clothing/eye protection/face protection.
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plasticware.
- Store all reagents at 2–8°C. Bring the reagents to 18–26°C prior to use, and return to 2–8°C following use.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional and national regulations.
- Care should be taken to prevent the contamination of kit components.
- Do not use components past their expiration dates and do not mix components from different serials.
- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing and washing throughout this procedure are necessary to maintain precision and accuracy.
- The two controls must be used for each test series.
- Use only distilled or deionized water for preparation of the reagents used in the test.
- Unused microtiter wells should be stored sealed within the enclosed plastic bag at 2–8°C.
- For veterinary use only.

Preparation of Samples

Dilute test samples five hundred fold (1:500) with Sample Diluent prior to being assayed (e.g., by diluting 1 μl of sample with 500 μl of Sample Diluent).

Note: DO NOT DILUTE CONTROLS. Be sure to change tips for each sample. Samples must be thoroughly mixed prior to dispensing into the Coated Plate.

Test Procedure

All reagents must be allowed to come to 18–26°C before use then mixed by inverting and swirling. Use a separate pipette tip for each sample

1. Obtain antigen-coated plate(s) and record the sample position on a worksheet.
2. Dispense 100 μl of undiluted Negative Control into two appropriate wells.
3. Dispense 100 μl of undiluted Positive Control into two appropriate wells.
4. Dispense 100 μl of diluted sample into applicable wells.
5. Incubate for 30 minutes (± 2 min.) at 18–26°C.
6. Wash each well with approximately 350 μl of distilled or deionized water 3–5 times.
7. Dispense 100 μl of Conjugate into each well.
8. Incubate for 30 minutes (± 2 min.) at 18–26°C.
9. Repeat step 6.
10. Dispense 100 μl of TMB Substrate solution into each well.
11. Incubate for 15 minutes (± 1 min.) at 18–26°C.
12. Dispense 100 μl of Stop Solution into each well to stop the reaction.
13. Blank the reader with air.
14. Measure and record absorbance values at 650 nm, A(650).

Results

For the assay to be valid, the difference between the Positive Control mean and the Negative Control mean ($PC\bar{x} - NC\bar{x}$) should be greater than 0.075. The Negative Control mean absorbance should be less than or equal to 0.150. The presence or absence of antibody to APV is determined by relating the A(650) value of the sample to the Positive Control mean. The Positive Control is standardized and represents significant antibody levels to APV in chicken and turkey serum. The relative level of antibody in the sample is determined by calculating the sample to positive (S/P) ratio. Endpoint titers are calculated using the equation described in the calculations section.

NOTE: IDEXX has instrument and software systems available which calculate means and % values and provide data summaries.

Calculation

Calculation of Negative Control mean

(NC \bar{x})

$$\text{NC}\bar{x} = \frac{\text{NC1 A650} + \text{NC2 A650}}{2}$$

Example:

$$\text{NC}\bar{x} = \frac{0.105 + 0.095}{2} = 0.100$$

Calculation of Positive Control mean

(PC \bar{x})

$$\text{PC}\bar{x} = \frac{\text{PC1 A650} + \text{PC2 A650}}{2}$$

Example:

$$\text{PC}\bar{x} = \frac{0.295 + 0.305}{2} = 0.300$$

S/P calculation for each sample

$$\text{S/P} = \frac{\text{Sample Mean} - \text{NC}\bar{x}}{\text{PC}\bar{x} - \text{NC}\bar{x}}$$

Example:

$$\text{S/P} = \frac{0.425 - 0.100}{0.300 - 0.100} = 0.162$$

Titer-Relates S/P at a 1:500 dilution to an endpoint titer:

$$\text{Log}_{10}\text{Titer} = 1.09 (\text{Log}_{10} \text{S/P}) + 3.36$$

Interpretation of Results

- Serum samples with S/P ratios of less than or equal to 0.20 should be considered negative.
- Serum samples with S/P ratios greater than 0.20 (titers greater than 396) should be considered positive and indicate vaccination or other exposure to APV.

Each laboratory should establish its own criterion for immunity with respect to antibody titer based on correlation of IDEXX APV to current laboratory test methodologies and on historical antibody responses to specific vaccines and vaccination protocols.

The immune status of a flock is best assessed by monitoring and recording antibody titers in representative samples as a function of time.

The resulting flock profiles allow an assessment of the distribution of antibody titers and an analysis of changes in titer over time.

👁 = *Modification in the using instructions*

Summarized Test Procedure

IDEXX strongly recommends that you read the complete instructions carefully before using the test the first time.

Steps	Action
1. Preparation and distribution of samples	Dilute test samples five hundred fold (1:500) with Sample Diluent prior to being assayed (e.g., by diluting 1 μ l of sample with 500 μ l of Sample Diluent). Note: DO NOT DILUTE CONTROLS. Be sure to change tips for each sample. Samples must be thoroughly mixed prior to dispensing into the coated plate. Dispense 100 μ l of undiluted Negative Control into two appropriate wells. Dispense 100 μ l of undiluted Positive Control into two appropriate wells. Dispense 100 μ l of diluted sample into applicable wells.
2. Sample incubation	Incubate for 30 minutes (\pm 2 min.) at 18–26°C.
3. Washing the plate	Wash each well with approximately 350 μ l of distilled or deionized water 3–5 times.
4. Conjugate Distribution	Dispense 100 μ l of Conjugate into each well.
5. Conjugate Incubation	Incubate for 30 minutes (\pm 2 min.) at 18–26°C .
6. Repeat step 4	
7. Substrate distribution	Dispense 100 μ l of TMB Substrate into each well.
8. Substrate incubation	Incubate for 15 minutes (\pm 1 min.) at 18–26°C.
9. Stopping the reaction	Dispense 100 μ l of Stop Solution into each well to stop the reaction.
10. Measure the plate	Measure and record absorbance values at 650 nm, A(650).
11. Interpretation	Serum samples with S/P ratios of less than or equal to 0.20 should be considered negative. Serum samples with S/P ratios greater than 0.20 (titers greater than 396) should be considered positive and indicate vaccination or other exposure to APV.

Manufactured by:

IDEXX Switzerland AG
Stationsstrasse 12
3097 Liebefeld-Bern, Switzerland

For technical assistance:

Contact your IDEXX area manager or distributor
or visit: www.idexx.com/production/contact
IDEXX Technical Support: 00-800-727-43399

*IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

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IDEXX Europe B.V.