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Chicken Anemia Virus Antibody Test Kit

For veterinary use only.

Name and Intended Use

IDEXX CAV is IDEXX's enzyme immunoassay for the detection of antibody to Chicken Anemia Virus (CAV) in chicken serum.

General Information

Chicken Anemia Virus (CAV) is an important pathogen of poultry and has been found in broilers, breeders and SPF flocks. Virus isolation is difficult and time consuming. Screening for the presence of antibody will indicate exposure to the virus. The enzyme-linked immunosorbent assay has been utilized to detect antibody against CAV. This method is quite useful in large-scale testing of flocks for exposure to CAV. The IDEXX CAV assay can be used as a screen for the presence of Chicken Anemia Virus antibodies or to assess response to vaccination.

Description and Principles

The IDEXX CAV assay is performed in a microtiter well coated with Chicken Anemia Virus antigen. During the first incubation, CAV antibodies present in the sample react with immobilized antigens. After a wash step, an Anti-CAV monoclonal antibody enzyme conjugate is added to the microwell. If no Anti-CAV antibodies are present in the test sample, the Anti-CAV conjugate is free to react with the CAV antigen. Conversely, if Anti-CAV antibodies are present in the sample, the enzyme-conjugated monoclonal antibodies are blocked from reacting with the antigen. Following this incubation period, the unreacted conjugate is removed by washing and a substrate/chromogen solution is added. In the presence of enzyme, the substrate is converted to a product which reacts with the chromophore to generate a blue color. The absorbance at 650nm, A(650), is measured using a spectrophotometer. Results are calculated by dividing the A(650) of the sample by the mean A(650) of the negative control, resulting in an S/N value. The quantity of antibodies to CAV is inversely proportional to the A(650) and thus, to the S/N value. The presence of CAV antibodies indicates previous exposure to Chicken Anemia Virus.

Reagent		Volume
1	CAV Antigen Coated Plates; preserved with Kathon	5
2	Positive Control — diluted chicken anti-CAV serum; preserved with sodium azide	2 mL
3	Negative Control — diluted chicken serum non-reactive to CAV; preserved with sodium azide	2 mL
4	Conjugate — anti-CAV: HRPO conjugate; preserved with gentamicin and Kathon	50 mL
5	Sample Diluent — preserved with sodium azide	235 mL
A	TMB Substrate	60 mL

B	Stop Solution	60 mL
C	Wash Concentrate (10X) — preserved with gentamicin	235 mL
Other Components: Zip lock bag.		1

NOTE: see table on page 21 for the description of international symbols used on the kit labels.

Materials Required but Not Provided

Precision pipettes and multiple delivery pipetting device with disposable pipette tips, 96-well plate reader, tubes for diluting samples, distilled or deionized water and device for the delivery and aspiration of wash solution. Reagent volumes listed in the Test Procedure require pipette precision of $\pm 5\%$.

Precautions and Warnings for Users

Handle all CAV biological materials as though capable of transmitting CAV. The antigen coated plates may be a source of CAV. Prior to coating on the solid phase, the antigen has been inactivated by chemical treatment. Nevertheless, do not assume complete inactivation. Some kit components contain sodium azide as a preservative. Dispose of contents in accordance with local, regional, and national regulations. Do not expose TMB solutions to strong light or any oxidation agents. Store all reagents at 2–8°C. All wastes should be properly decontaminated prior to disposal. Do not use kit serials past expiration date and do not intermix components from kits with different serial numbers. Careful pipetting and washing throughout this procedure are necessary to maintain precision and accuracy. Optimal results will be obtained by strict adherence to this protocol. For veterinary use only.

Preparation of Samples

Note: Dilution of serum samples may be made with microdilution tubes or directly in the microtiter plate by adding the appropriate amount of serum to the previously pipetted Sample Diluent. Be sure to change tips for each sample. Samples must be thoroughly mixed prior to dispensing into the coated plate. **DO NOT DILUTE THE CONTROLS.**

A. Field Exposure

For detection of antibody subsequent to field exposure to Chicken Anemia Virus. Dilute serum samples ten-fold (1/10) with sample diluent prior to being assayed (e.g. diluting 10 μL of sample with 90 μL of sample diluent).

B. Vaccination

For the detection of elevated antibody levels (e.g. subsequent to CAV vaccination), dilute serum samples 1/100 (e.g. diluting 10 μL of sample with 990 μL of sample diluent). This dilution will allow adequate differentiation of responses from samples within a vaccinated flock. Responses of vaccination will depend upon dose, route of administration, age of flock, etc.

Preparation of Wash Solution

The Wash Concentrate (10X) should be brought to 18–26°C and mixed to ensure dissolution of any precipitated salts. The Wash Concentrate (10X) must be diluted 1/10 with distilled/deionized water before use (e.g. 30 mL of concentrate plus 270 mL of water per plate to be assayed).

Test Procedure

Allow the reagents to come to 18–26°C, then mix gently by inverting and swirling.

1. Obtain antigen coated plate(s) and record the sample position. If using partial plates, remove only those wells sufficient for samples to be tested. Place the remaining wells, along with the desiccant, in the extra ziplock bag provided and return to 2–8°C.
2. Dispense 100 μL of UNDILUTED Negative Control into duplicate wells.
3. Dispense 100 μL of UNDILUTED Positive Control into duplicate wells.
4. Dispense 100 μL of diluted sample into appropriate wells. Samples may be tested in duplicate but a single well is acceptable.
5. Incubate for 60 minutes (\pm 5 minutes) at 18–26°C.
6. Wash each well with approximately 350 μL of Wash Solution 3-5 times.
7. Dispense 100 μL of Conjugate into each well.
8. Incubate for 30 minutes (\pm 2 minutes) at 18–26°C.
9. Repeat step 6.
10. Dispense 100 μL of TMB Substrate solution into each well.
11. Incubate for 15 minutes (\pm 1 minute) at 18–26°C.
12. Dispense 100 μL of Stop Solution into each well to stop the reaction.
13. Measure and record absorbance values at 650 nm, A(650).

Results

For the assay to be valid, the Negative Control optical density (650 nm) must be greater than or equal to 0.600 and the Positive Control S/N must be less than or equal to 0.50. For invalid tests, technique may be suspect and the assay should be repeated. The presence or absence of antibody to CAV is determined by the sample to negative (S/N) ratio for each sample.

Calculations

1.	Negative Control mean ($\text{NC}\bar{x}$)	$\frac{\text{NC1 A}(650) + \text{NC2 A}(650)}{2} = \text{NC}\bar{x}$
2.	Positive Control mean ($\text{PC}\bar{x}$)	$\frac{\text{PC1 A}(650) + \text{PC2 A}(650)}{2} = \text{PC}\bar{x}$
3.	S/N Ratio	$\frac{\text{Sample A}(650)}{\text{NC}\bar{x}} = \text{S/N}$

Interpretation of Results

1/10 Dilution

1. Samples with S/N ratios greater than 0.60 are considered negative within the limits of the test.
2. Samples with S/N ratios less than or equal to 0.60 are considered positive.

1/100 Dilution

The immune status of a vaccinated flock is best assessed by monitoring and recording antibody levels in representative samples as a function of time. The resulting flock profiles allow an assessment of the distribution of antibody levels and an analysis of changes in status over time.

IDEXX Technical Services:

IDEXX USA Tel: 1 800 548 9997 or 1 207 556 4890

Fax: 1 800 328 5461 or 1 207 556 4826

IDEXX Europe Tel: 00800 727 43399

Fax: 00800 433 99329

U.S. Vet. License No. 313

Product Code: 50C1.00

*IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

©2011 IDEXX Laboratories, Inc. All rights reserved.

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© 2011 IDEXX Laboratories, Inc. Alle Rechte vorbehalten.
IDEXX Europe B.V.