**Classical Swine Fever Virus (CSFV)
Antibody Test Kit**

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Test With Confidence*

Classical Swine Fever Virus (CSFV) Antibody Test Kit

For veterinary use only

Name and Intended Use

IDEXX CSFV Ab is IDEXX's enzyme immunoassay for the detection of CSFV specific antibodies in swine serum or plasma.

General Information

Classical Swine Fever Virus (CSFV), Bovine Viral Diarrhoea Virus (BVDV), and Border Disease Virus (BDV) are the 3 members of the genus Pestivirus within the family Flaviviridae. CSFV causes serious losses in the pig industry since it is highly pathogenic and can cause widespread deaths. Pigs infected with highly virulent CSFV strains may shed a high amount of virus before showing clinical signs of the disease. Animals that survive an acute or sub-acute infection develop antibodies and will no longer spread the virus. Moderately virulent, less pathogenic strains may lead to chronic infection, when pigs excrete the virus continuously or intermittently until death. Fetuses of susceptible pregnant sows, exposed to CSFV, might become infected through the placenta. Congenital infection can result in abortion, fetal mummification, stillborn and/or weak piglets or embryonic malformations. The most frequent outcome of congenital infection with low virulent strains is the birth of persistently infected piglets, spreading the virus without signs of the disease, in the absence of immune response. Pigs may also be infected with BVDV or BDV. Although these infections in pigs are usually mild and self-limiting, it is important to reliably distinguish them from the CSFV infection.

Descriptions and Principles

IDEXX CSFV Ab is IDEXX's enzyme immunoassay for the detection of CSFV specific antibodies in swine serum or plasma. The assay is a blocking ELISA which utilizes microplates coated with CSFV antigen. Antibodies present in the test sample will block the binding of Horseradish Peroxidase conjugated CSFV specific monoclonal antibodies. The bound monoclonal antibodies are detected by a substrate reactive with Horseradish Peroxidase. The result is indicated by color development. The optical density is measured by a microplate reader at a single wavelength of 450 nm, or dual wavelengths of 450 nm and 650 nm. Color development is weak (positive result), when CSFV specific antibodies are present in the test sample. Color development is maximal (negative result) in the absence of specific antibodies. The blocking percentage of the sample is calculated from the optical density (OD450) obtained with the test sample and the OD450 of the negative control.

Reagents

Store all reagents at 2–8°C.

Reagents		Volume	Volume
1	CSFV Antigen Coated Plates	2	5
2	Positive Control	1 ml	1 ml
3	Negative Control	1 ml	1 ml
4	(Anti-CSFV HRPO) Conjugate	25 ml	60 ml
5	Sample Diluent	15 ml	45 ml
A	TMB Substrate N.12	20 ml	60 ml
B	Stop solution N.3	20 ml	60 ml
C	Wash concentrate (10X)	125 ml	480 ml

NOTE: see table on page 32 for the description of international symbols used on the kit labels.

Materials Required but Not Provided

- Precision pipettes or multiple delivery pipetting devices suitable for delivering 10 to 1000 μ l. Reagents volumes listed in the “Test Procedure” require pipette precision less than or equal to 5%.
- Disposable pipette tips
- 500 ml graduated cylinder for wash solution
- 96-well Microplate reader equipped with a single wavelength filter of 450 nm or a dual wavelength of 450 nm and 650 nm
- Ultrapure water (Milli Q or double distilled)
- Device for the delivery and aspiration of wash solution
- A trap for retaining aspirate and disinfectant
- Humid chamber and plate covers or plate sealers
- Vortex mixer
- Centrifuge tubes or micro tubes

Precautions and Warnings for Users

- Handle all biological material as potentially infectious.
- Do not pipette by mouth.
- Do not eat, drink or smoke where specimens or kit reagents are being handled.
- The substrate solution is irritating to the eyes, respiratory system and skin. Avoid contact with skin and eyes. Wear protective gloves/protective clothing/eye protection/face protection
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plasticware.
- Store all reagents at 2–8°C. Bring the reagents to 18–26°C prior to use, and return to 2–8°C following use.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional and national regulations
- Care should be taken to prevent the contamination of kit components.
- Do not use components past their expiration dates and do not mix components from different serials.
- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing and washing throughout this procedure are necessary to maintain precision and accuracy.
- The two controls must be used for each test series.
- Use only distilled or deionized water for preparation of the reagents used in the test.
- Unused microtiter wells should be stored sealed within the enclosed plastic bag at 2–8°C.
- For veterinary use only.

Preparation of Reagents

Coated microtiter plates, sample diluent, positive control, negative control, conjugate, substrate solution and stop solution are ready-to-use.

Wash Solution

The Wash Concentrate (10X) should be brought to 18–26°C and mixed to ensure dissolution of any precipitated salts. The Wash Concentrate (10X) must be diluted 1:10 in ultrapure water (Milli Q or double distilled). When prepared under sterile conditions, the reconstituted Wash Solution can be stored for one week at 2–8°C. Example: For 300 ml Wash Solution, dilute 30 ml Wash Concentrate (10X) with 270 ml ultrapure water (Milli Q or double distilled) and mix well.

Preparation of Samples

Fresh, refrigerated (maximum of 8 days at 2–8°C), or frozen serum or plasma can be tested.

Test Procedure

All reagents must be allowed to come to 18–26°C before use.

Reagents should be mixed by gentle swirling or vortexing. Use a separate pipette tip for each sample.

1. Add 50 μ l of the Sample Diluent to all wells to be tested and to the control wells.
2. Add 50 μ l of the Positive Control and Negative Control to the appropriate wells. Use a separate pipette tip for different control sera.
3. Add 50 μ l of the serum or plasma sample to the remainder of the plate. Use a separate pipette tip for each sample.
4. Mix the contents of the microwells by gently tapping the plate or use a shaker for microtiter plates.
5. Incubate for 2 hours (\pm 5 min.) or overnight (12-18 hours) at 18–26°C. With either option, the plates should be tightly sealed or incubated in a humid chamber using plate covers to avoid any evaporation.
6. Wash each well with approximately 300 μ l of Wash Solution three times. Aspirate liquid contents of all wells after each wash. Following the final aspiration, firmly tap residual wash fluid from each plate onto absorbent material. Avoid plate drying between washes and prior to the addition of the next reagent.
7. Add 100 μ l of conjugate to each well. Incubate 30 minutes (\pm 2 min.) at 18–26°C sealed or in a humid chamber using plate covers.
8. Repeat Step 6.
9. Add 100 μ l of TMB Substrate N.12 to each well and incubate 10 minutes (\pm 1 min.) at 18–26°C away from direct light. Begin timing after the first well is filled.
10. Stop the reaction at 10 minutes by adding 100 μ l of Stop Solution N.3 to each well. Add the stop solution in the same order as the substrate solution was added in step 9.
11. Measure the absorbance of the samples and controls at 450 nm, or using a dual wavelength of 450 nm and 650 nm on a microplate reader (use air as blank).
12. Calculate the mean absorbance values (see 'Calculations') for each test sample and for the controls.

Results

For the assay to be valid, the Negative Control mean OD450 should be greater than 0.500. The blocking percentage of the Positive Control should be greater than 50%. For invalid assays, technique may be suspect and the assay should be repeated following a thorough review of the package insert.

NOTE: IDEXX has instrument and software systems available which calculate means and % blocking and provide data summaries.

Calculation

Calculation of Negative

Control mean (NC \bar{x})

$$NC\bar{x} = \frac{NC1 A_{450} + NC2 A_{450}}{2}$$

Example:

$$NC\bar{x} = \frac{1.280 + 1.300}{2} = 1.290$$

Calculation of Positive

Control mean (PC \bar{x})

$$PC\bar{x} = \frac{PC1 A_{450} + PC2 A_{450}}{2}$$

Example:

$$PC\bar{x} = \frac{0.120 + 0.140}{2} = 0.130$$

Calculation for test samples

Blocking-percentage:

$$\frac{NC\bar{x} A_{450} - \text{Sample } A_{450}}{NC\bar{x} A_{450}} \times 100$$

Example: sample result 1.219

$$\frac{1.290 - 1.219}{1.290} \times 100 = 5.5\%$$

Interpretation of Results

- The test sample is positive (antibodies are present) if it gives a blocking percentage $\geq 40\%$.
- The test sample is negative (antibodies are absent) if it gives a blocking percentage $\leq 30\%$.
- If the blocking percentage of the sample is $> 30\%$ and $< 40\%$, retest the animal at a later date. If the sample is still doubtful after retesting, the result can be confirmed by Serum Neutralization Test. Please refer to the regulations currently in force in your country.

☞ = *Modification in the using instructions*

Summarized Test Procedure

IDEXX strongly recommends that you read the complete instructions carefully before using the test the first time.

Steps	Action							
1. Sample diluent distribution	All reagents should equilibrate to 18– 26°C before use. Reagents should be mixed by gentle swirling or vortexing. Add 50 μ l of the Sample Diluent to all wells to be tested and to the control wells.							
2. Sample distribution	Add 50 μ l of the Positive Control and Negative Control to the appropriate wells. Use a separate pipette tip for different control sera. Add 50 μ l of the serum or plasma sample to the remainder of the plate. Use a separate pipette tip for each sample. Mix the contents of the microwells by gently tapping the plate or use a shaker for Microtiter plates.							
3. Sample incubation	Incubate for 2 hours (\pm 5 min.) or overnight (12-18 hours) at 18–26°C. With either option, the plates should be tightly sealed or incubated in a humid chamber using plate covers to avoid any evaporation.							
4. Washing the plate	Aspirate liquid contents of all wells into an appropriate waste reservoir. Wash each well with approximately 300 μ l of Wash Solution three times. Aspirate liquid contents of all wells after each wash. Following the final aspiration, firmly tap residual wash fluid from each plate onto absorbent material.							
5. Conjugate Distribution	Add 100 μ l of Conjugate to each well.							
6. Conjugate Incubation	Incubate 30 (\pm 2 min.) minutes at 18– 26°C sealed or in a humid chamber using plate covers.							
7. Repeat step 4								
8. Substrate distribution	Add 100 μ l of TMB Substrate N.12 into each well.							
9. Substrate incubation	Incubate 10 minutes (\pm 1 min.) at 18– 26°C away from direct light. Begin timing after the first well is filled.							
10. Stopping the reaction	Stop the reaction at 10 minutes by adding 100 μ l of Stop Solution N.3 to each well. Add the stop solution in the same order as the substrate solution was added in step 8.							
11. Measure the plate	Measure the absorbance of the samples and controls at 450 nm, or using a dual wavelength of 450 nm and 650 nm on a microplate reader (use air as blank). Calculate the results.							
12. Interpretation %	<table border="1"> <tbody> <tr> <td rowspan="2">Serum, Plasma</td> <td>\leq30%</td> <td>> 30% to <40%</td> <td>\geq40%</td> </tr> <tr> <td>Negative</td> <td>Suspect</td> <td>Positive</td> </tr> </tbody> </table>	Serum, Plasma	\leq 30%	> 30% to <40%	\geq 40%	Negative	Suspect	Positive
Serum, Plasma	\leq 30%		> 30% to <40%	\geq 40%				
	Negative	Suspect	Positive					

Manufactured by:

IDEXX Switzerland AG
Stationsstrasse 12
3097 Liebefeld-Bern, Switzerland

For technical assistance:

Contact your IDEXX area manager or distributor
or visit: www.idexx.com/production/contact
IDEXX Technical Support: 00-800-727-43399

*IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

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IDEXX Europe B.V.