**Classical Swine Fever Virus Antigen Test Kit
(CSFV) / Serum Plus**

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Classical Swine Fever Virus Antigen Test Kit (CSFV) / Serum Plus

For veterinary use only

Name and Intended Use

IDEXX CSFV Ag Serum Plus is IDEXX's enzyme immunoassay for the detection of CSFV antigens in swine serum, plasma and tissue (preferably tonsils, spleen, mesenteric lymph node or kidneys) samples.

General Information

Classical Swine Fever Virus (CSFV), Bovine Viral Diarrhoea Virus (BVDV), and Border Disease Virus (BDV) are the 3 members of the genus Pestivirus within the family *Flaviviridae*. CSFV causes serious losses in the pig industry since it is an obligate pathogen and can cause widespread deaths. Pigs infected with highly virulent CSFV strains may shed a high amount of virus before showing clinical signs of the disease. Animals that survive an acute or sub-acute infection develop antibodies and will no longer spread the virus. Moderately virulent, less pathogenic strains may lead to chronic infection, when pigs excrete the virus continuously or intermittently until death. Fetuses of susceptible pregnant sows, exposed to CSFV, might become infected through the placenta. Congenital infection can result in abortion, fetal mummification, stillborn and/or weak piglets or embryonic malformations. The most frequent outcome of congenital infection with low virulent strains is the birth of persistently infected piglets, spreading the virus without signs of the disease, in the absence of immune response. Pigs may also be infected with BVDV or BDV. Although these infections in pigs are usually mild and self limiting, virus detection in such cases is very rare. Since the IDEXX CSFV Ag Serum Plus test has high sensitivity for all pestiviruses, the detection of BVDV or BDV cannot be excluded. Positive results should therefore be confirmed with a CSFV specific assay.

Descriptions and Principles

A microtitration format has been configured by immobilizing specific monoclonal antibodies on the plates. CSFV Ag of the sample is captured on the plates. After incubation of the test sample in the well, captured CSFV Ag is detected by specific polyclonal antiserum for CSFV (E^{ms}) and a horseradish-peroxidase Conjugate. Next, unbound Conjugate is washed away and a Substrate/chromogen solution is added. In the presence of enzyme, Substrate is converted into a product which reacts with the chromogen to generate a blue color. Upon addition of the Stop Solution, a yellow color is generated. The absorbance at a single wavelength of 450 nm [A(450)] or a dual wavelength of 450 nm and 650 nm [A(450/650)] is measured using a spectrophotometer. The corrected OD value of the samples is calculated by using the absorbance [A(450)] or [A(450/650)] obtained with the test sample and corrected for the absorbance of the Negative Control.

Reagents

		Volume	
1	Anti-E ^{ms} mAb Coated Plate	2	5
2	Positive Control	1 x 1.6 mL	1 x 2.0 mL
3	Negative Control	1 x 1.6 mL	1 x 2.0 mL
4	Conjugate	1 x 25 mL	1 x 60 mL
5	Tissue Soaking Buffer Concentrate (2X)	1 x 115 mL	2 x 115 mL

A	TMB Substrate N.12	1 x 20 mL	1 x 60 mL
B	Stop Solution N.3	1 x 20 mL	1 x 60 mL
C	Wash Concentrate (10X)	1 x 125 mL	1 x 480 mL
D	Detection Solution	1 x 15 mL	1 x 30 mL

Note: See table at the end of the insert for a description of international symbols used on the labels of this kit.

Storage

Store the reagents at 2–8°C. Reagents are stable until expiration date, provided they have been stored properly.

Materials Required but Not Provided

- Precision micropipettes or multi-dispensing micropipettes
- Disposable pipette tips
- Graduated cylinder for wash solution
- 96-well microplate reader (equipped with single 450 nm filter or dual filters at 450 nm and 650 nm)
- Microplate washer (manual, semi-automatic or automatic system)
- Use only distilled or deionized water for preparation of the reagents used in the test
- Microplate covers (lid, aluminium foil or adhesive)
- Humid Chamber/ Incubator capable of maintaining a temperature of +37°C ($\pm 3^{\circ}\text{C}$)
- Vortex or equivalent
- Centrifuge (capacity 1500 x g)

Precautions and Warnings for Users

- Handle all biological material as potentially infectious.
- The Substrate solution is irritating to eyes, respiratory system, and skin. Avoid contact with skin and eyes.
- Wear protective gloves/protective clothing/eye or face protection when handling samples and reagents.
- Refer to the product Material Safety Data Sheet for additional information.

Laboratory Practices

- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing, and washing throughout this procedure are necessary to maintain precision and accuracy.
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plastic ware.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional, and national regulations.
- Care should be taken to prevent contamination of kit components. Do not pour unused reagents back into containers.
- Do not use kit past expiration date and do not intermix components from kits with different lot numbers.

Preparation of Reagents

Wash Solution

The Wash Concentrate (10X) should be brought to 18–26°C and mixed to ensure dissolution of any precipitated salts. The Wash Concentrate (10X) must be diluted 1:10 with distilled/deionized water before use (e.g., 30 mL of Wash Concentrate (10X) plus 270 mL of distilled/deionized water per plate to be assayed). When prepared under sterile conditions, the Wash Solution can be stored for one week at 2–8°C.

Tissue soaking buffer

The Tissue Soaking Buffer Concentrate (2X) should be brought to 18–26°C and mixed to ensure dissolution of any precipitated salts. The Tissue Soaking Buffer Concentrate (2X) must be diluted 1:2 with distilled/deionized water before use. The amount of diluted Tissue Soaking Buffer needs to be calculated as required. When prepared under sterile conditions, the diluted Tissue Soaking Buffer can be stored for up to two weeks at 2–8°C.

Sample preparation

Serum and Plasma Samples

Fresh or frozen serum and plasma samples can be tested.

Tissues

Use preferably fresh tissues but, if necessary, tissues can be stored frozen. Process one or two tissues from each animal submitted, preferably tonsils, spleen, mesenteric lymph node or kidneys.

- a.) With scissors, cut tissue into small pieces (250 mg).
- b.) Place the tissue pieces into an appropriate centrifuge tube or an Eppendorf tube and add 1 mL Tissue Soaking Buffer. Vortex. Leave at 18–26°C for 1–2 hours.
- c.) Centrifuge at 1,500 x g for 10 minutes and use 50 μ L of the clean supernatant for testing as described in the Test Procedure.

Please use the overnight incubation protocol for tissue samples!

Test Procedure

All reagents must be allowed to come to 18–26°C before use. Reagents should be mixed by gentle swirling or vortexing. Use a separate pipette tip for each sample.

1. Obtain coated plate(s) and record the sample position on a worksheet.
2. Add 50 μ L of Detection Solution into each well. A multichannel pipette (8 or 12 channel) can be used for this step.
3. Add 50 μ L of Negative Control into two appropriate wells.
4. Add 50 μ L of Positive Control into two appropriate wells.
5. Add 50 μ L of samples into remaining wells.
6. Mix the content of the microwells by gently tapping the plate or use a shaker for microtiter plates.
7. Incubate samples:
 - a) Serum and plasma samples: incubate for 2 hours (\pm 5 min.) at +37°C (\pm 3°C) or 12–18 hours at 2–8°C. With either option, plates should be tightly sealed or incubated in a humid chamber using plate covers to avoid any evaporation.
 - b) Tissue samples: incubate overnight (12–18 hours) at 2–8°C. The plates should be tightly sealed or incubated in a humid chamber using plate covers to avoid any evaporation.

- 8. Remove the solution and wash each well with approximately 300 μL of Wash Solution 5 times. Avoid plate drying between plate washings and prior to the addition of the next reagent. Tap each plate onto absorbent material after the final wash to remove any residual wash fluid.
Important! Control carefully that no traces of blood are left on the walls or edges of the wells. Additional 2–3 washes can be necessary to remove the blood before proceeding to the next step.
- 9. Dispense 100 μL of Conjugate into each well.
- 10. Incubate for 30 minutes (± 3 min.) at 18–26°C.
- 11. Repeat step 8.
- 12. Dispense 100 μL of TMB Substrate N.12 solution into each well.
- 13. Incubate 10 minutes (± 1 min.) at 18–26°C away from direct light. Begin timing after the first well is filled.
- 14. Dispense 100 μL of Stop Solution N.3 into each well to stop the reaction. Add the Stop Solution in the same order as the Substrate solution.
- 15. Blank the spectrophotometer on air.
- 16. Measure and record the absorbance of the samples and controls at 450 nm or using a dual wavelength of 450 nm and 650 nm.
- 17. Calculate the results.

Results

For the assay to be valid, the difference (P–N) between the Positive Control mean (PC \bar{x}) and the Negative Control mean (NC \bar{x}) must be greater than or equal to 0.150 optical density (OD). In addition, the Negative Control mean (NC \bar{x}) must be less than or equal to 0.250 OD. For invalid assays, technique may be suspect and the assay should be repeated following a thorough review of the package insert. The presence or absence of CSFV antigen in the sample is determined by corrected OD value (S–N) for each sample.

See CALCULATIONS for examples.

Note: IDEXX has instrument and software systems available that calculate means and % values and provide data summaries.

Calculation

Negative Control mean (NC \bar{x})	$\text{NC}\bar{x} = \frac{\text{NC1 A450} + \text{NC2 A450}}{2}$
Positive Control mean (PC \bar{x})	$\text{PC}\bar{x} = \frac{\text{PC1 A450} + \text{PC2 A450}}{2}$
S–N for test samples	$\text{S–N} = \text{Sample A450} - \text{NC}\bar{x}$

Interpretation of Results

- Samples with S–N values equal or less than 0.300 are classified as negative for CSFV Ag.
- Samples with S–N values greater than 0.300 are considered positive for CSFV Ag.

☞ = Major change in the user instructions

Summarized Test Procedure

IDEXX strongly recommends that you read the complete instructions carefully before using the test the first time.

Steps	Action
1. Preparation of reagents	The Wash Concentrate (10X) must be diluted 1:10 with distilled/deionized water before use. The Tissue Soaking Buffer Concentrate (2X) must be diluted 1:2 with distilled/deionized water before use.
2. Sample distribution (Serum, plasma or tissue samples)	Add 50 μ L of Detection Solution into each well. A multichannel pipette (8 or 12 channel) can be used for this step. Add 50 μ L of Negative Control into two appropriate wells. Add 50 μ L of Positive Control into two appropriate wells. Add 50 μ L of samples into remaining wells. Mix the content of the microwells by gently tapping the plate or use a shaker for microtiter plates.
3. Sample Incubation	a. Serum and plasma samples: Incubate samples for 2 hours (\pm 5 min.) at 37°C (\pm 3°C) or overnight (12–18 hours) at 2–8°C. With either option, the plates should be tightly sealed or incubated in a humid chamber using plate covers to avoid any evaporation. b. Tissue samples: Incubate samples overnight (12–18 hours) at 2–8°C. The plates should be tightly sealed or incubated in a humid chamber using plate covers to avoid any evaporation.
4. Washing the plate	Remove the solution and wash each well with approximately 300 μ L of Wash Solution 5 times.
5. Conjugate distribution	Dispense 100 μ L of Conjugate into each well.
6. Conjugate incubation	Incubate for 30 minutes (\pm 3 min.) at 18–26°C.
7. Repeat step 4	
8. Substrate distribution	Dispense 100 μ L of TMB Substrate N.12 solution into each well.
9. Substrate incubation	Incubate 10 minutes (\pm 1 min.) at 18–26°C away from direct light. Begin timing after the first well is filled.
10. Stopping the reaction	Dispense 100 μ L of Stop Solution N.3 into each well to stop the reaction.
11. Measure the plate	Blank the spectrophotometer on air. Measure and record the absorbance of the samples and controls at 450 nm or using a dual wavelength of 450 nm and 650 nm.
12. Interpretation (S–N) (Serum, plasma and tissue samples)	Samples with S–N values equal or less than 0.300 are classified as negative for CSFV Ag. Samples with S–N values greater than 0.300 are considered positive for CSFV Ag.

For technical assistance:

Contact your IDEXX area manager or distributor or visit our website: www.idexx.com/production/contact
IDEXX Technical Support: +800 727 43399

*IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

©2013 IDEXX Laboratories, Inc. All rights reserved.

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©2013 IDEXX Laboratories, Inc. Alle Rechte vorbehalten.
IDEXX Switzerland AG
IDEXX Europe B.V.