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Test With Confidence*

Influenza A Virus Antibody Test Kit

For veterinary use only

Name and Intended Use

IDEXX Influenza A Test Kit is an enzyme immunoassay for the detection of antibody to influenza A virus in animal serum.

General Information

Domestic and wild animal species are affected by the influenza A virus. The disease is characterized by a wide range of responses that include virtually no clinical signs to variable mortality. The severity of influenza A will vary with the strain and/or isolate of the virus, the age and the immune status of the animal, the presence of concurrent viral infections and whether influenza A is complicated by secondary infections. Respiratory signs include coughing, nasal and/or ocular discharge, sneezing and dyspnea. Hyperthermia may be observed with associated anorexia, weight loss, lethargy and prostration. Because of the variation and severity of clinical symptoms, serological testing has significant advantages. Monitoring for exposure of a flock or herd to influenza A is facilitated by the measurement of antibody to influenza A in serum.

Descriptions and Principles

This assay is designed to measure the relative level of antibody to influenza A in animal serum. The assay is performed on 96-well plates that have been coated with influenza A viral antigen. Upon incubation of the test sample in the coated wells, influenza A specific antibody forms a complex with the coated antigen. After washing away unbound material, an anti-influenza A nucleoprotein monoclonal antibody enzyme conjugate is added to the wells. In the absence of influenza A antibodies in the test sample, the conjugate is free to bind to the Influenza A antigen on the plate. Conversely, if there are antibodies to influenza A present in the sample, the anti-influenza A nucleoprotein conjugate is blocked from binding to the antigen. Unbound conjugate is washed away and enzyme substrate is added. Subsequent color development is inversely proportional to the amount of anti-influenza A antibodies in the test sample.

Reagents

Store all reagents at 2–8°C.

Reagents		Volume
1	Influenza A Antigen Coated Plate	5
2	Positive Control	2.0 mL
3	Negative Control	2.0 mL
4	Conjugate (Anti-Influenza A: HRPO)	55 mL
5	Dilution Buffer	235 mL
A	TMB Substrate	60 mL
B	Stop Solution	60 mL
C	Wash Concentrate (10X)	235 mL

NOTE: see table on page 28 for the description of international symbols used on the kit labels.

Materials Required but Not Provided

- Precision pipettes or multiple delivery pipetting devices (reagents volumes listed in the “Test Protocol” require pipette precision of $\leq 5\%$)
- Disposable pipette tips
- 96-well plate reader equipped with 650 nm filter
- Tubes for diluting samples
- Distilled or deionized water
- Device for the delivery and aspiration of wash solution

Precautions and Warnings for Users

- Handle all influenza A biological materials as though capable of transmitting influenza A. The antigen-coated plates may be a source of influenza A. Prior to coating on the solid phase, the antigen has been inactivated by chemical treatment. Nevertheless, do not assume complete inactivation.
- Some kit components contain sodium azide as a preservative. Disposal requires flushing plumbing with large volumes of water to prevent formation of copper or lead azide complexes, which may explode upon percussion. All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional and national regulations.
- Do not expose TMB solutions to strong light or any oxidation agents. Store all reagents at 2–8°C
- Do not use components past expiration date and do not intermix components from kits with different lot numbers.
- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing and washing throughout this procedure are necessary to maintain precision and accuracy.
- For veterinary use only.
- Wear protective gloves/protective clothing/eye protection/face protection.

Preparation of Reagents

Wash Solution

The Wash Concentrate (10X) should be brought to 18–26°C and mixed to ensure dissolution of any precipitated salts. The Wash Concentrate (10X) must be diluted 1 to 10 with distilled/deionized water before use (e.g., 30 ml of Wash Concentrate (10X) plus 270 ml of water per plate to be assayed). When prepared under sterile conditions, the Wash Solution can be stored for one week at 2–8°C.

Preparation of Samples

Dilute test samples tenfold (1/10) with the Dilution Buffer prior to being assayed (e.g. by diluting 15 μL of sample with 135 μL of Dilution Buffer).

DO NOT DILUTE THE CONTROLS.

Use a separate pipette tips for each sample. Samples should be mixed thoroughly before being distributed on the plates.

Test Procedure

All reagents must be allowed to come to 18–26°C before use.

Reagents should be mixed by gentle swirling or vortexing. Use a separate pipette tip for each sample.

1. Obtain coated plates and record the sample position on a worksheet.
2. Dispense 100 μL of UNDILUTED Negative Control into two wells.
3. Dispense 100 μL of UNDILUTED Positive Control into two wells.
4. Dispense 100 μL of diluted samples into remaining wells.
5. Incubate for 60 minutes (± 5 min.) at 18–26°C.
6. Wash each well with approximately 350 μL of Wash Solution 3–5 times.
7. Dispense 100 μL of Conjugate into each well.
8. Incubate for 30 minutes (± 2 min.) at 18–26°C.
9. Repeat step 6.
10. Dispense 100 μL of TMB Substrate solution into each well.
11. Incubate 15 minutes (± 1 min.) at 18–26°C.
12. Dispense 100 μL of Stop Solution into each well to stop the reaction.
13. Blank the spectrophotometer on air.
14. Measure and record the absorbance of the samples and controls at 650 nm.
15. Calculate the results.

Results

For the assay to be valid, the Negative Control mean ($\text{NC}\bar{x}$) must be greater than or equal to 0.600 optical density (OD). In addition, the Positive Control mean S/N must be less than 0.50. For invalid assays, technique may be suspect and the assay should be repeated following a thorough review of the package insert. The presence or absence of antibody to Influenza A is determined by the sample to negative (S/N) ratio for each sample.

Note: IDEXX has instrument and software systems available which calculate means and %blocking and provide data summaries.

Calculation

Calculation of Negative

Control mean (NC \bar{x})

$$NC\bar{x} = \frac{NC1 \text{ A650} + NC2 \text{ A650}}{2}$$

Calculation of Positive Control

mean (PC \bar{x})

$$PC\bar{x} = \frac{PC1 \text{ A650} + PC2 \text{ A650}}{2}$$

Calculation for test sample

$$S/N = \frac{\text{Sample mean A650}}{NC\bar{x}}$$

Interpretation of Results

Avian Cutoff:

- Samples with S/N values < 0.50 should be considered influenza A antibody positive.
- Samples with an S/N ratio \geq 0.50 are considered negative for the presence of influenza A antibodies.

Other Animal Species Cutoff:

- Samples with S/N values < 0.60 should be considered influenza A antibody positive.
- Samples with an S/N ratio \geq 0.60 are considered negative for the presence of influenza A antibodies.

Note: Heat inactivation of sera from geese may give a weak false positive test result.

👁 = *Modification in the using instructions*

Summarized Test Procedure

IDEXX strongly recommends that you read the complete instructions carefully before using the test the first time.

Steps	Action
1. Preparation of Reagents	The Wash Concentrate (10X) must be diluted 1 to 10 with distilled/deionized water before use.
2. Sample Preparation	Dilute test samples tenfold (1/10) with the Dilution Buffer prior to being assayed.
3. Sample Distribution	Dispense 100 μ L of UNDILUTED Negative Control into two wells. Dispense 100 μ L of UNDILUTED Positive Control into two wells. Dispense 100 μ L of diluted samples into remaining wells.
4. Sample incubation	Incubate for 60 minutes (\pm 5 min.) at 18–26°C.
5. Washing the plate	Wash each well with approximately 350 μ L of Wash Solution 3–5 times.
6. Conjugate distribution	Dispense 100 μ L of Conjugate into each well.
7. Conjugate incubation	Incubate for 30 minutes (\pm 2 min) at 18–26°C
8. Repeat step 5	
9. Substrate Distribution	Dispense 100 μ L of TMB Substrate solution into each well
10. Substrate Incubation	Incubate 15 minutes (\pm 1 min) at 18–26°C
11. Stopping the reaction	Dispense 100 μ L of Stop Solution into each well to stop the reaction
12. Measure the plate	Blank the spectrophotometer on air. Measure and record the absorbance of the samples and controls at 650 nm. Calculate the results.
13. Interpretation	<u>Avian Cutoff:</u> Samples with S/N values < 0.50 should be considered influenza A antibody positive. Samples with an S/N ratio \geq 0.50 are considered negative for the presence of influenza A antibodies. <u>Other Animal Species Cutoff:</u> Samples with S/N values < 0.60 should be considered influenza A antibody positive. Samples with an S/N ratio \geq 0.60 are considered negative for the presence of influenza A antibodies.

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For technical assistance:

Contact your IDEXX area manager or distributor
or visit: www.idexx.com/production/contact
IDEXX Technical Support: 00-800-727-43399

*IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

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IDEXX Montpellier SAS
IDEXX Europe B.V.