**Mycoplasma Gallisepticum Antibody Test Kit**

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For veterinary use only.

### Name and Intended Use

IDEXX MG is IDEXX's enzyme immunoassay for the detection of antibody to *Mycoplasma gallisepticum* (Mg) in chicken and turkey serum.

### General Information

An assessment of immune status as well as serologic identification of Mg requires a measurement of antibody to Mg in serum. Enzyme immunoassay systems have proven efficacious in the quantification of antibody levels to Mg and facilitate the monitoring of immune status in large flocks.

### Descriptions and Principles

This assay is designed to measure the relative level of antibody to Mg in chicken and turkey serum. Mg antigen is coated on 96-well plates. Upon incubation of the test sample in the coated well, antibody specific to Mg forms a complex with the coated antigens. After washing away unbound material from the wells, a conjugate is added which binds to any attached antibody in the wells. Unbound conjugate is washed away and enzyme substrate is added. Subsequent color development is directly related to the amount of antibody to Mg present in the test sample.

Reagent		Volume
1	Mg Antigen Coated Plates	5
2	Positive Control — diluted chicken anti-Mg, preserved with sodium azide	1 x 1.9 mL
3	Negative Control — diluted chicken serum non-reactive to Mg, preserved with sodium azide	1 x 1.9 mL
4	Conjugate — (Goat) anti-chicken/anti-turkey: HRPO Conjugate; preserved with gentamicin and Kathon	1 x 50 mL
5	Sample Diluent — buffer, preserved with sodium azide	1 x 235 mL
A	TMB Substrate	1 x 60 mL
B	Stop Solution	1 x 60 mL

**NOTE:** see table on page 21 for the description of international symbols used on the kit labels.

## Materials Required but Not Provided

Precision pipettes and multiple delivery pipetting device with disposable pipette tips, 96-well plate reader, tubes for diluting samples, distilled or deionized water and device for the delivery and aspiration of wash solution. Reagent volumes listed in the Test Procedure require pipette precision of  $\pm 5\%$ .

## Precautions and Warnings for Users

Handle all Mg biological materials as though capable of transmitting Mg. The antigen coated plates may be a source of Mg. Prior to coating on the solid phase, the antigen has been inactivated by chemical treatment. Nevertheless, do not assume complete inactivation. Some kit components contain sodium azide as a preservative. Dispose of contents in accordance with local, regional, and national regulations. Do not expose TMB solutions to strong light or any oxidation agents. Store all reagents at 2–8°C. All wastes should be properly decontaminated prior to disposal. Do not use kit serials past expiration date and do not intermix components from kits with different serial numbers. Careful pipetting and washing throughout this procedure are necessary to maintain precision and accuracy. Optimal results will be obtained by strict adherence to this protocol. For veterinary use only.

## Preparation of Samples

Dilute test samples five hundred fold (1:500) with sample diluent prior to being assayed (e.g., by diluting 1  $\mu\text{L}$  of sample with 500  $\mu\text{L}$  of Sample Diluent). **NOTE: DO NOT DILUTE CONTROLS. Be sure to change tips for each sample. Samples must be thoroughly mixed prior to dispensing into the coated plate.**

## Test Procedure

Allow the reagents to come to 18–26°C, then mix gently by inverting and swirling.

1. Obtain antigen-coated plate(s) and record the sample.
2. Dispense 100  $\mu\text{L}$  of UNDILUTED Negative Control into duplicate wells.
3. Dispense 100  $\mu\text{L}$  of UNDILUTED Positive Control into duplicate wells.
4. Dispense 100  $\mu\text{L}$  of diluted sample into appropriate wells. Samples may be tested in duplicate but a single well is acceptable.
5. Incubate for 30 minutes ( $\pm 2$  minutes) at 18–26°C.
6. Wash each well with approximately 350  $\mu\text{L}$  of distilled or deionized water 3–5 times. Aspirate completely.
7. Dispense 100  $\mu\text{L}$  of Conjugate into each well.
8. Incubate for 30 ( $\pm 2$  minutes) minutes at 18–26°C.
9. Repeat step 6.
10. Dispense 100  $\mu\text{L}$  of TMB Substrate Solution into each well.
11. Incubate for 15 minutes ( $\pm 1$  minute) at 18–26°C.
12. Dispense 100  $\mu\text{L}$  of Stop Solution into each well to stop the reaction.
13. Measure and record absorbance values at 650nm, A(650).

## Results

For the assay to be valid, the difference between the Positive Control mean and the Negative Control mean ( $PC\bar{x} - NC\bar{x}$ ) should be greater than 0.075. The Negative Control mean absorbance should be less than or equal to 0.150. The presence or absence of antibody to Mg is determined by relating the A(650) value of the unknown to the Positive Control mean. The Positive Control is standardized and represents significant antibody levels to Mg in serum. The relative level of antibody in the sample is determined by calculating the sample to positive (S/P) ratio. Endpoint titers are calculated using the equation described in the calculations section.

## Calculations

1.	Negative Control mean ( $NC\bar{x}$ )	$\frac{NC1 A(650) + NC2 A(650)}{2} = NC\bar{x}$
2.	Positive Control mean ( $PC\bar{x}$ )	$\frac{PC1 A(650) + PC2 A(650)}{2} = PC\bar{x}$
3.	S/P Ratio	$\frac{\text{Sample Mean} - NC\bar{x}}{PC\bar{x} - NC\bar{x}} = S/P$
4	Titer - Relates S/P at a 1:500 dilution to an endpoint titer:	$\text{Log}_{10}\text{Titer} = 1.09 (\log_{10} S/P) + 3.36$

## Interpretation of Results

Serum samples with S/P ratios of less than or equal to 0.50 should be considered negative. S/P ratios greater than 0.50 (titers greater than 1076) should be considered positive and indicates exposure to Mg. Each laboratory should establish its own criterion for immunity with respect to antibody titer based on correlation of IDEXX MG to current laboratory test methodologies and on historical antibody responses.

**IDEXX Technical Services**

IDEXX USA Tel: 1 800 548 9997 or 1 207 556 4890

Fax: 1 800 328 5461 or 1 207 556 4826

IDEXX Europe Tel: 00800 727 43399

Fax: 00800 433 99329

U.S. Vet. License No. 313  
Product Code 5070.00

\*IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

©2011 IDEXX Laboratories, Inc. All rights reserved.

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© 2011 IDEXX Laboratories, Inc. Alle Rechte vorbehalten.
IDEXX Europe B.V.