## **Newcastle Disease Virus Antibody Test Kit**

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## Newcastle Disease Virus Antibody Test Kit

(For use in Chicken)

For veterinary use only.

### Name and Intended Use

IDEXX NDV is IDEXX's enzyme immunoassay for the detection of antibody to Newcastle disease virus (NDV) in chicken serum.

### General Information

An assessment of immune status as well as serologic identification of NDV requires a measurement of antibody to NDV in serum. Enzyme immunoassay systems have proven efficacious in the quantification of antibody levels to NDV and facilitate the monitoring of immune status in large flocks.

### Descriptions and Principles

This assay is designed to measure the relative level of antibody to NDV in chicken serum. Viral antigen is coated on 96-well plates. Upon incubation of the test sample in the coated well, antibody specific to NDV forms a complex with the coated viral antigens. After washing away unbound material from the wells, a conjugate is added which binds to any attached chicken antibody in the wells. Unbound conjugate is washed away and enzyme substrate is added. Subsequent color development is directly related to the amount of antibody to NDV present in the test sample.

Reagents	Volume	
1	NDV Antigen Coated Plate	5
2	Positive Control — diluted chicken anti-NDV serum; preserved with sodium azide	1 x 1.9 mL
3	Negative Control — diluted chicken serum non-reactive to NDV; preserved with sodium azide	1 x 1.9 mL
4	Conjugate — (Goat) anti-chicken: HRPO conjugate; preserved with gentamicin and Kathon	1 x 50 mL
5	Sample Diluent — buffer, preserved with sodium azide	1 x 235 mL
A	TMB Substrate	1 x 60 mL
B	Stop Solution	1 x 60 mL

**NOTE:** See table at the end of the insert for a description of international symbols used on the labels of this kit.

## Materials Required but Not Provided

Precision pipettes and multiple delivery pipetting device with disposable pipette tips, 96-well plate reader, tubes for diluting samples, distilled or deionized water and device for the delivery and aspiration of wash solution. Reagent volumes listed in the Test Procedure require pipette precision less than or equal to 5%.

## Precautions and Warnings for Users

Handle all NDV biological materials as though capable of transmitting NDV. The antigen-coated plates may be a source of NDV. Prior to coating on the solid phase, the antigen has been inactivated by chemical treatment. Nevertheless, do not assume complete inactivation. Some kit components contain sodium azide as a preservative. Dispose of contents in accordance with local, regional, and national regulations. Do not expose TMB solutions to strong light or any oxidation agents. Store all reagents at 2–8°C. All wastes should be properly decontaminated prior to disposal. Do not use kit serials past expiration date and do not intermix components from kits with different serial numbers. Careful pipetting and washing throughout this procedure are necessary to maintain precision and accuracy. Optimal results will be obtained by strict adherence to this protocol. For veterinary use only.

## Preparation of Samples

Dilute test samples five hundred fold (1:500) with sample diluent prior to being assayed (e.g., by diluting 1  $\mu\text{L}$  of sample with 500  $\mu\text{L}$  of Sample Diluent). **NOTE: DO NOT DILUTE CONTROLS.** Be sure to change tips for each sample. Samples must be thoroughly mixed prior to dispensing into the coated plate.

## Test Procedure

Allow the reagents to come to 18–26°C, then mix gently by inverting and swirling.

1. Obtain antigen-coated plate(s) and record the sample position.
2. Dispense 100  $\mu\text{L}$  of UNDILUTED Negative Control into duplicate wells.
3. Dispense 100  $\mu\text{L}$  of UNDILUTED Positive Control into duplicate wells.
4. Dispense 100  $\mu\text{L}$  of diluted sample into appropriate wells. Samples may be tested in duplicate but a single well is acceptable.
5. Incubate for 30 minutes ( $\pm$  2 minutes) at 18–26°C.
6. Aspirate liquid content of all wells into appropriate waste reservoir.
7. Wash each well with approximately 350  $\mu\text{L}$  of distilled or deionized water 3–5 times. Aspirate completely.
8. Dispense 100  $\mu\text{L}$  of Conjugate into each well.
9. Incubate for 30 minutes ( $\pm$  2 minutes) at 18–26°C.
10. Repeat steps 6 and 7.
11. Dispense 100  $\mu\text{L}$  of TMB substrate into each well.
12. Incubate for 15 minutes ( $\pm$  1 minute) at 18–26°C.
13. Dispense 100  $\mu\text{L}$  of Stop Solution into each well to stop the reaction.
14. Measure and record absorbance values at 650nm, A(650).

## Results

For the assay to be valid, the difference between the Positive Control mean and the Negative Control mean ( $PC\bar{x} - NC\bar{x}$ ) should be greater than 0.075. The Negative Control mean absorbance should be less than or equal to 0.150. The presence or absence of antibody to NDV is determined by relating the A(650) value of the unknown to the Positive Control mean. The Positive Control is standardized and represents significant antibody levels to NDV in chicken serum. The relative level of antibody in the sample is determined by calculating the sample to positive (S/P) ratio. Endpoint titers are calculated using the equation described in the calculations section.

## Calculations

1	Negative Control mean ( $NC\bar{x}$ )	$\frac{NC1 A(650) + NC2 A(650)}{2} = NC\bar{x}$
2	Positive Control mean ( $PC\bar{x}$ )	$\frac{PC1 A(650) + PC2 A(650)}{2} = PC\bar{x}$
3	S/P Ratio	$\frac{\text{Sample Mean} - NC\bar{x}}{PC\bar{x} - NC\bar{x}} = S/P$
4	Titer - Relates S/P at a 1:500 dilution to an endpoint titer:	$\text{Log}_{10} \text{Titer} = 1.09 (\log_{10} S/P) + 3.36$

## Interpretation of Results

Serum samples with S/P ratios of less than or equal to 0.20 should be considered negative. S/P ratios greater than 0.20 (titers greater than 396) should be considered positive and indicate vaccination or other exposure to NDV. Each laboratory should establish its own criterion for immunity with respect to antibody titer based on correlation of IDEXX NDV to current laboratory test methodologies and on historical antibody responses to specific vaccines and vaccination protocols. The immune status of a flock is best assessed by monitoring and recording antibody titers in representative samples as a function of time. The resulting flock profiles allow an assessment of the distribution of antibody titers and an analysis of changes in titer over time.

## **IDEXX Technical Services**

IDEXX USA Tel: 1 800 548 9997 or 1 207 556 4890

Fax: 1 800 328 5461 or 1 207 556 4826

IDEXX Europe Tel: 00800 727 43399

Fax: 00800 433 99329

U.S. Vet. License No. 313

Product Code: 5080.00

\*IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

©2012 IDEXX Laboratories, Inc. All rights reserved.

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© 2012 IDEXX Laboratories, Inc. Alle Rechte vorbehalten.
IDEXX Europe B.V.