## **Porcine Reproductive and Respiratory Syndrome Virus Antibody Test Kit for Oral Fluids**

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## Porcine Reproductive and Respiratory Syndrome Virus Antibody Test Kit for Oral fluids

For veterinary use only.

### Name and Intended Use

IDEXX PRRS OF is an enzyme immunoassay for the detection of antibodies to porcine reproductive and respiratory syndrome virus (PRRSV) in samples of swine oral fluid collected from individual pigs or composite samples from a pen of pigs\*.

### General Information

Porcine reproductive & respiratory syndrome virus (PRRSV) is an arterivirus that was first reported in both Europe and the US in the late 1980's as distinct strains with variable antigenicity<sup>1</sup>. The virus causes reproductive problems (abortions, weak piglets, failure to farrow, and sow deaths), respiratory disease (pneumonia, difficulty breathing, poor performance, and death) and mild neurologic signs. PRRS is considered to be one of the most economically significant diseases affecting the swine industry. An assessment of exposure to the PRRSV as a result of natural infection or vaccination is facilitated by a measurement of antibodies in serum or plasma. Recently, detection of antibodies in oral fluid samples of pigs is increasingly being used as a tool to monitor PRRS exposure in pigs<sup>2,3</sup>. It is important that the immunodiagnostic assay measuring the antibodies be sensitive to the different strains of the virus. For veterinary use only.

### Descriptions and Principles

IDEXX PRRS OF is IDEXX's enzyme immunoassay for the detection antibody to porcine reproductive and respiratory syndrome virus (PRRSV) in samples of swine oral fluid. A microtitration format has been configured by coating recombinant PRRSV antigens on the plate. Upon incubation of the test sample in the coated well, antibody specific to PRRSV forms a complex with the coated antigens. After washing away unbound material from the wells, an anti-porcine:horseradish peroxidase conjugate is added that binds to any porcine antibody attached in the wells. Unbound conjugate is washed away and TMB substrate is added to the wells. Color development is then stopped using an acid stop solution, which turns the color from blue to yellow. The intensity of the color is proportional to the amount of bound specific antibodies against PRRSV present in the sample.

\*Claim for performance on pen samples is based on a study of 500 pigs divided into pens of 25 pigs each. Pigs vaccinated with modified live vaccine (Ingelvac® PRRS MLV, Boehringer Ingelheim Vetmedica Inc.) 14 days prior to the day of the study (at peak of sero-conversion for antibodies) were introduced into pens of pigs sero-negative for PRRS antibodies at different prevalence of antibody-positive pigs (0%, 4%, 12%, 20%, and 36%), with 5 pens for each of the prevalence values. All pigs were bled for confirmation of serum antibody titers. Oral fluids were collected from each pen by hanging a cotton rope on the pen for 30 minutes; five successive 30-minute replicate collections were done per pen by hanging a separate new rope per collection event. PRRS antibodies were detected in at least 96% of all collection events at a prevalence of at least 20% sero-positive (vaccinated) pigs.

Reagents		Volume
1	PRRSV Antigen Coated Plate	5
2	Positive Control	1 x 4 mL
3	Negative Control	1 x 4 mL
4	Conjugate	1 x 60 mL
5	Sample Diluent	1 x 120 mL
A	TMB Substrate N.12	1 x 60 mL
B	Stop Solution N.3	1 x 60 mL
C	Wash Concentrate (10X)	1 x 235 mL
<b>Other Components:</b> Zip lock bag.		1

**Note:** see table at the end of the insert for a description of international symbols used on the labels of this kit.

### Materials Required but Not Provided

- Precision pipettes or multiple-delivery pipetting devices suitable for delivering 5 to 1000  $\mu$ L (reagents volumes listed in the “Test Procedure” require pipette precision less than or equal to 5%).
- Microplate covers (lid, aluminium foil or adhesive)
- Disposable pipette tips
- Graduated cylinder for wash solution
- 96-well Microplate reader equipped with a dual 450nm and 650nm filter
- Distilled or deionized water
- Device for the delivery and aspiration of wash solution
- A trap for retaining aspirate and disinfectant
- Tubes for diluting samples

### Precautions and Warnings for Users

- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing and washing throughout this procedure are necessary to maintain precision and accuracy.
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plastic ware.
- Store all reagents at 2–8°C. Bring the reagents to 18–26°C prior to use, and return to 2–8°C following use.
- Care should be taken to prevent contamination of kit components.
- Do not use kits past their expiration dates and do not intermix components from kits with different serials.
- Use only distilled or deionized water for preparation of the reagents used in the test.
- Unused microtiter wells should be stored sealed in a plastic bag with desiccant at 2–8°C.
- The two controls must be used for each test series.

## Sample Collection and Handling

Samples of oral fluids can be collected from individual pigs or collectively from pens using hanging ropes according to widely publicized methods available online or in the literature<sup>2,3</sup>.

## Warning

When testing oral fluids from pigs on transition diets in the first weeks after weaning, exercise caution upon interpreting results, as oral fluids from pigs feeding on diets supplemented with sprayed-dried plasma of porcine origin may test positive due to the presence of PRRSV antibodies in the porcine plasma supplement<sup>4</sup>.

## Laboratory Practices

- Do not pipette by mouth.
- Do not eat, drink, or smoke where specimens or kit reagents are being handled.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional and national regulations.
- The Substrate solution is irritating to eyes, respiratory system, and skin. Avoid contact with skin and eyes.

## Preparation of Reagents

### Wash Solution

Determine the amount of Wash Solution needed for washing the microtiter plates. Dilute the Wash Concentrate (10X) 1:10 with distilled/deionized water (1 part concentrate with 9 parts water, e.g., 30 mL Wash concentrate (10X) + 270 mL distilled water). When prepared under sterile conditions (sterile water and sterile flasks or containers) this solution can be used during seven days when stored between 2–8°C.

### Preparation of samples

Dilute samples to be tested 1:2 with the Sample Diluent (e.g., by diluting 100  $\mu$ L of sample with 100  $\mu$ L of Sample Diluent). **Note:** Do not dilute controls. Be sure to change tips for each sample and record the position of each sample on the plate using a worksheet. Samples should be mixed prior to dispensing into the PRRSV coated wells.

## Test Procedure

All reagents must be allowed to come to 18–26°C before use then mixed by inverting and swirling. Use a separate pipette tip for each sample.

1. Obtain antigen-coated plate(s) and record the sample position.
2. Dispense 100  $\mu$ L of UNDILUTED Negative Control into two wells of the assay plate.
3. Dispense 100  $\mu$ L of UNDILUTED Positive Control into two wells of the assay plate.
4. Dispense 100  $\mu$ L of diluted samples into appropriate wells.
5. Cover the wells and incubate for 2 hours ( $\pm$ 5 min.) at 18–26°C.
6. Empty or aspirate liquid content from the microwells and wash each well with approximately 300  $\mu$ L of Wash Solution three to five times. Discard liquid contents of all wells after each wash. Avoid plate drying between plate washings and prior to the addition of conjugate. Following the final fluid removal, gently, but firmly tap residual wash fluid from each plate onto absorbent material.

7. Dispense 100  $\mu$ L of Conjugate into each well.
8. Incubate for 30 minutes ( $\pm 2$  min.) at 18–26°C.
9. Repeat step 6.
10. Dispense 100  $\mu$ L of TMB Substrate N.12 into each test plate well.
11. Incubate for 15 minutes ( $\pm 1$  min.) at 18–26°C.
12. Dispense 100  $\mu$ L of Stop Solution N.3 into each well of the test plate to stop the reaction.  
Add the Stop Solution N.3 in the same order as the Substrate solution was added in step 10.
13. Measure and record the A(450 nm-REF) for samples and controls within 30 minutes of dispensing the Stop solution N.3. REF is the reference wavelength, which can be within the range of A(620 nm) to A(650 nm).
14. Calculate results.

## Results

For the assay to be valid, the following specifications must be met: the Positive Control mean minus the mean of the Negative Control ( $PC\bar{x} - NC\bar{x}$ ) must be greater than or equal to 0.150. In addition, the Negative Control ( $NC\bar{x}$ ) must be less than or equal to 0.150. For invalid assays, technique may be suspect and the assay should be repeated following a thorough review of the product insert. The presence or absence of antibody to PRRSV is determined by calculating the sample to positive (S/P) ratio. See Calculation section for examples.

**Note:** IDEXX has instrument and software systems available that calculate means and S/P and provide data summaries.

## Calculations

Examples:

Negative Control mean ( $NC\bar{x}$ )	$NC\bar{x} = \frac{NC1 A(450-REF) + NC2 A(450-REF)}{2}$	$\frac{0.064 + 0.08}{2} = 0.072$
Positive Control mean ( $PC\bar{x}$ )	$PC\bar{x} = \frac{PC1 A(450-REF) + PC2 A(450-REF)}{2}$	$\frac{0.560 + 0.596}{2} = 0.578$
Calculation for test sample	$S/P = \frac{Sample - NC\bar{x}}{PC\bar{x} - NC\bar{x}}$	Sample = A(450-REF) = 0.750 $S/P = \frac{0.750 - 0.072}{0.578 - 0.072} = 1.34$

## Interpretation of Results

The presence or absence of antibody to PRRSV is determined by calculating the S/P ratio for each sample.

- If the S/P ratio is less than 0.40, the sample is classified as NEGATIVE for PRRSV antibodies.
- If the S/P ratio is greater than or equal to 0.40, then the sample is classified as POSITIVE for PRRSV antibodies.

## Summarized Test Procedure

IDEXX strongly recommends that you read the complete instructions carefully before using the test the first time.

Step	Action
<b>1. Preparation of Reagents</b>	Determine the amount of Wash Solution needed for washing the microtiter plates. Dilute the Wash Concentrate (10X) 1:10 with distilled/deionized water (1 part concentrate with 9 parts water, e.g., 30 mL Wash Concentrate (10X) + 270 mL distilled water).
<b>2. Preparation of Samples</b>	Dilute samples to be tested 1:2 with the Sample Diluent (e.g., by diluting 100 $\mu$ L of sample with 100 $\mu$ L of Sample Diluent).
<b>3. Sample Distribution</b>	Dispense 100 $\mu$ L of UNDILUTED Negative Control into two wells of the assay plate. Dispense 100 $\mu$ L of UNDILUTED Positive Control into two wells of the assay plate. Dispense 100 $\mu$ L of diluted samples into appropriate wells.
<b>4. Sample Incubation</b>	Cover the wells and incubate for 2 hours ( $\pm$ 5 min.) at 18–26°C.
<b>5. Washing the Plate</b>	Empty or aspirate liquid contents of all wells and wash each well with approximately 300 $\mu$ L of Wash Solution three to five times. Discard liquid contents of all wells after each wash. Avoid plate drying between plate washings and prior to the addition of conjugate. Following the final liquid removal, gently, but firmly tap residual wash fluid from each plate onto absorbent material.
<b>6. Conjugate Distribution</b>	Dispense 100 $\mu$ L of Conjugate into each well.
<b>7. Conjugate Incubation</b>	Incubate for 30 minutes ( $\pm$ 2 min.) at 18–26°C.
<b>8. Repeat step 5</b>	
<b>9. Substrate Distribution</b>	Dispense 100 $\mu$ L of TMB Substrate N.12 into each test plate well.
<b>10. Substrate Incubation</b>	Incubate for 15 minutes ( $\pm$ 1 min.) at 18–26°C.
<b>11. Stop the reaction</b>	Dispense 100 $\mu$ L of Stop Solution N.3 into each well of the test plate to stop the reaction.
<b>12. Measure the plate</b>	Measure and record the A(450-REF) for samples and controls within 30 minutes of stopping the reaction. Calculate results.
<b>13. Interpretation</b>	If the S/P ratio is less than 0.40, the sample is classified as NEGATIVE for PRRSV antibodies. If the S/P ratio is greater than or equal to 0.40, then the sample is classified as POSITIVE for PRRSV antibodies.

## References

1. J. Zimmerman and K.-J. Yoon, Editors, The PRRS Compendium (second ed.), National Pork Board, Des Moines, IA (2003).
2. Kittawornrat A, et al. Detection of Porcine reproductive and respiratory syndrome virus (PRRSV) antibodies in oral fluid specimens using a commercial PRRSV serum antibody enzyme-linked immunosorbent assay. *J Vet Diagn Invest*, 24(2):262-9 (2012).
3. Prickett JR and Zimmerman JJ. The development of oral fluid-based diagnostics and applications in veterinary medicine. *An Health Res Rev*, 11: 207-216 (2010).
4. Johnson JK, Main R, Zimmerman J. 2012. Exogenous sources of PRRSV antibody detectable in the PRRSV oral fluid ELISA. *J Swine Health Prod* (in press).

## For technical assistance:

Contact your IDEXX area manager or distributor or visit our website:[www.idexx.com/production/contact](http://www.idexx.com/production/contact)  
IDEXX Technical Support: 00 800 727 43399

\*IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

©2013 IDEXX Laboratories, Inc. All rights reserved.

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© 2013 IDEXX Laboratories, Inc. Alle Rechte vorbehalten.
IDEXX Switzerland AG
IDEXX Europe B.V.