## **Pseudorabies Virus gpl Antibody Test Kit**

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Test With Confidence™

**IDEXX**

## Pseudorabies Virus gpl Antibody Test Kit

For veterinary use only.

### Name and Intended Use

IDEXX PRV/ADV gI is an enzyme immunoassay for the detection of antibodies in swine serum to the gpl antigen of pseudorabies virus (PRV/Aujeszky's Disease). The presence of antibodies to gpl indicates exposure to field strains of PRV and/or vaccines containing gpl antigen. In the United States, it is intended for use in management and pseudorabies eradication applications. When used with gpl-deleted PRV vaccines, for example Intervet- Schering Plough Animal Health, and when performed at approved laboratories, it is an official differential pseudorabies test for use in the Cooperative State/Federal Pseudorabies Eradication Program.

### General Information

Strategies traditionally used to provide protection for swine from the effects of pseudorabies infection have included the administration of conventional and gene-deleted, USDA-licensed PRV vaccines.<sup>1</sup> Standard serologic procedures used to assess PRV exposure following vaccination are of limited value as they cannot distinguish naturally infected from vaccinated swine. Safe and efficacious gpl-deleted PRV vaccines that allow for serologic differentiation have been developed by the manufacturers listed above. The virus used in these vaccines was selected such that the synthesis of virulence factors were eliminated without compromising the immunogenicity of the virus. In addition, a mechanism for serologic differentiation is provided by the natural deletion of the DNA sequences that code for gpl, an apparently non-essential viral protein. Thus, the IDEXX PRV/ADV gI assay, being specific for antibodies to gpl, ignores antibody titers in animals vaccinated with gpl-deleted PRV vaccines by approved manufacturers, but detects animals infected with field strains or vaccinated with strains that contain the gpl antigen. The assay uses monoclonal antibodies specific for gpl. Prior to vaccination with gpl-deleted products, swine should be tested using either this assay or the IDEXX PRV/ADV gB assays to determine immune status. Subsequent to vaccination, the swine should be routinely monitored at least semiannually for exposure to field strains of PRV using the IDEXX PRV/ADV gI assay.

## Descriptions and Principles

The IDEXX PRV/ADV gI assay is performed in a PRV antigen-coated microwell using a two-fold (1:2) serum dilution. During the first incubation, PRV antibodies present in the serum, including those produced against gpl, react with antigens on the plate. Subsequent to a wash step, an anti-PRV-gpl monoclonal antibody conjugate is added to the microwell and is allowed to compete for the gpl viral antigen during a second incubation. If no gpl antibodies are present in the test serum, the conjugated gpl antibodies are free to react with the gpl antigen. Conversely, if gpl antibodies are present in the test serum, the enzyme-conjugated monoclonal antibodies are blocked from reacting with the antigen. Following this incubation period, the unreacted conjugate is removed by washing, and a substrate/chromogen solution is added. In the presence of enzyme, the substrate is converted to a product that reacts with the chromophore to generate a blue color. The absorbance at 650 nm, A(650), is measured using a spectrophotometer. Results are calculated by dividing the A(650) of the sample by the mean A(650) of the negative control, resulting in a sample/negative (S/N) value. The quantity of antibodies to gpl is inversely proportional to the A(650) and, thus, to the S/N value. The presence of PRV antibodies, including anti-gpl, indicates a previous exposure to a field strain of PRV or application of conventional modified live or killed virus vaccines. The presence of PRV antibodies detected by the IDEXX PRV/ADV gB assay, but absence of antibodies to gpl antigen as assessed by the IDEXX PRV/ADV gI assay, indicates a response to a gpl-deleted vaccine.

Reagents		Volume	
1	PRV Antigen Coated Plate	6	30
2	Positive Control — diluted porcine anti-PRV gpl serum; preserved with sodium azide	1 x 5.0 mL	1 x 5.0 mL
3	Negative Control — porcine serum non-reactive to PRV gpl; preserved with sodium azide	1 x 5.0 mL	1 x 5.0 mL
4	Conjugate — anti-PRV gpl; HRPO conjugate; preserved with gentamicin and Kathon	1 x 60 mL	1 x 350 mL
5	Sample Diluent — preserved with sodium azide	1 x 120 mL	1 x 300 mL
A	TMB Substrate	1 x 60 mL	1 x 315 mL
B	Stop Solution	1 x 60 mL	1 x 315 mL
C	Wash Concentrate (10X) — preserved with gentamicin	1 x 235 mL	3 x 480 mL
<b>Other Components:</b> Zip lock bag.		1	1

**Note:** See table at the end of the insert for a description of symbols used on the insert and labels of this kit.

## Storage

Store the reagents at 2–8°C. Reagents are stable until expiration date, provided they have been stored properly.

## Materials Required but Not Provided

- Precision micropipettes and multi-dispensing micropipettes
- Disposable pipette tips
- Graduated cylinder for wash solution
- 96-well microplate reader (equipped with 650 nm filter)
- Microplate washer (manual, semi-automatic or automatic system)
- Use only distilled or deionized water for preparation of the reagents used in the test
- Vortex or equivalent

## Precautions and Warnings

- Handle all biological material as potentially infectious. The antigen used in the kit reagents may not be completely inactivated.
- Wear protective gloves / protective clothing / eye or face protection when handling samples and reagents.
- Refer to the product Material Safety Data Sheet for additional information.
- See the end of this insert for reagent hazard and precaution warnings.

## Laboratory Practices

- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing, and washing throughout this procedure are necessary to maintain precision and accuracy. Use a separate pipette tip for each sample and control.
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plastic ware.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional, and national regulations.
- Care should be taken to prevent contamination of kit components. Do not pour unused reagents back into containers.
- Do not use kit past expiration date and do not intermix components from kits with different serial numbers.

## Preparation of Wash Solution

The Wash Concentrate (10X) should be brought to 18–26°C and mixed to assure dissolution of any precipitated salts. The Wash Concentrate (10X) must be diluted 1/10 with distilled/deionized water before use (e.g., 30 mL of concentrate plus 270 mL of water per plate to be assayed).

## Preparation of Samples

Dilute test samples two-fold (1:2) with Sample Diluent (e.g., by diluting 125  $\mu\text{L}$  of sample with 125  $\mu\text{L}$  of Sample Diluent). Be sure to change tips for each sample and record the position of each sample. Samples should be mixed prior to dispensing into the PRV-coated plate.

## Test Procedure

All reagents must be allowed to come to 18–26°C before use. Mix reagents by gentle inverting or swirling.

---

1 Obtain antigen-coated plate(s) and record the sample position. If using partial plates, remove only those wells sufficient for samples to be tested. Place the remaining wells, along with the desiccant, in the extra ziplock bag provided and return to 2–8°C.

---

2 Dispense 100  $\mu\text{L}$  of DILUTED Negative Control (NC) (DILUTED 1:2) into duplicate or triplicate wells.

---

3 Dispense 100  $\mu\text{L}$  of DILUTED Positive Control (PC) (DILUTED 1:2) into duplicate wells.

---

4 Dispense 100  $\mu\text{L}$  of DILUTED sample into appropriate wells.

---

5 Incubate for 60 minutes ( $\pm 5$  minutes) at 18–26°C or overnight at 2–8°C.

---

6 Remove the solution and wash each well with approximately 300  $\mu\text{L}$  of Wash Solution 3–5 times. Avoid plate drying between plate washings and prior to the addition of the next reagent. Tap each plate onto absorbent material after the final wash to remove any residual wash fluid.

---

7 Dispense 100  $\mu\text{L}$  of Conjugate into each well.

---

8 Incubate for 20 minutes ( $\pm 1$  minute) at 18–26°C.

---

9 Repeat Step 6.

---

10 Dispense 100  $\mu\text{L}$  of TMB Substrate into each well.

---

11 Incubate for 15 minutes ( $\pm 1$  minute) at 18–26°C.

---

12 Dispense 50  $\mu\text{L}$  of Stop Solution into each well.

---

13 Measure and record the A(650) for samples and controls.

## 14 Calculations:

### Controls

$$NC\bar{x} = \frac{A1 A(650) + A2 A(650) + A3 A(650)^*}{3}$$

$$PC\bar{x} = \frac{A4 A(650) + A5 A(650)}{2}$$

\*Example shows Negative Control run in Triplicate.

---

### Validity criteria

$$NC\bar{x} - PC\bar{x} \geq 0.300$$

For invalid assays, technique may be suspect and the assay should be repeated following a thorough review of the package insert.

---

### Samples

$$S/N = \frac{\text{Sample A}(650)}{NC\bar{x}}$$

The presence or absence of antibody to the antigen is determined by first calculating the S/N value of each sample.

---

## 15 Interpretation:

Negative	Suspect	Positive*
$S/N > 0.70$	$0.60 < S/N \leq 0.70$	$S/N \leq 0.60$

\*Confirm all positives in duplicate.

**Note:** IDEXX has instrument and software systems available which calculate results and provide data summaries.

## **Bibliography**

1. Gustafson DP, "Pseudorabies." Leman AD, Glock RD, Mengeling WL, Penny RH, Scholl E, Straw B, eds. Disease of Swine, 5th ed. Ames, Iowa: The Iowa State University Press;1984:209-223.

### **For technical assistance:**

IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4895

IDEXX Europe Tel: +800 727 43399

Contact your IDEXX area manager or distributor or visit our website: [idexx.com/contactlpd](http://idexx.com/contactlpd)

U.S. Vet. License No. 313

Product Code 5113.00

IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries

©2014 IDEXX Laboratories, Inc. All rights reserved.

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## WARNING / ATTENTION / ATENCIÓN / ATENÇÃO / ACHTUNG

### H316 / P332+P313 / EUH208

**Conjugate** – Causes mild skin irritation. If skin irritation occurs: Get medical advice/attention. Contains Kathon. May produce an allergic reaction.

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### H319 / P280 / P337+P313

**Sample Diluent** – Causes serious eye irritation. Wear eye protection/face protection. If eye irritation persists: Get medical advice/attention.

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**Stop solution** – Causes mild skin irritation. Causes serious eye irritation. Wear protective gloves/eye protection/face protection. If skin irritation occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention.

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IDEXX Europe B.V.